

# 6RUNNER® 10DLX (858WM)

SHOPRIDER®



Introducing a new member to Shoprider's 6Runner® series power chairs, the 6Runner® 10" Deluxe offers the enhanced stability of all 6 points of contact with the ground and true mid-wheel drive technology. The 6Runner® 10" Deluxe beautifully captures the functionality and maneuverability for someone on the go. Ride along with confidence with the exceptional quality you can always find in a Shoprider brand power chair.

## Available Colors



Red



Blue



## 6RUNNER® 10DLX (858WM)

- Luxurious leatherette seat
- Front and Rear suspension
- Stable four post seating system
- Reclining seat and adjustable seat height
- Adjustable length, height and angle foot plate
- Adjustable width, height and flip-up arm rest
- Midwheel drive maneuverability
- 10" wheels
- 300 lb. capacity



Adjustable height



Many seating options available



Flip-up arm rest



Reclining seat

## 6RUNNER® 10DLX SPECIFICATIONS

(858WM)

Overall Dimensions (L* <sup>1</sup> x W x H)	in	36x 22 x 38.6
Number, Size of Tires (front, middle, rear)	ea, in	2, 6"; 2, 10"; 2, 5"
Max. Single User Weight (on the level road)	lb	300
Total Unit Weight (w/battery, 50Ah* <sup>2</sup> )	lb	205
Battery Module Capacity	_V_Ah x Pcs	12V36 Ah up to 12V50Ah x 2
Charger	Amp	4A Off board
Output Power of Motor (for Reference only)	hp x Pcs	0.8 x 2
Turning Radius	in	18
Max. Speed* <sup>2</sup>	mph	5
Range* <sup>4</sup> (Per Charge w/STD 36Ah battery/300lbs. User Weight) ( After the battery and mechanical moving parts fully break in)	mile	20

All specification subject to change without prior notice.

**WARNING!** Please make sure that you are physically fit to operate this unit. Consult your dealer /doctor for advice.

\*1. Include the Anti-Tip Wheel.

\*2. Driver weight may exceed weight of the unit, speed must be reduced when turning.

\*3. When proceeding up any incline, please move your seat to the most forward position or if you have the deluxe seat, make sure that it is in the 90 degree (upright) position.

\*4. The actual driving range varies with the factors shown as below.

(A) The Weight of Occupant (E) Ambient Temperature  
(B) Level Ground (F) The Way of Driving  
(C) Battery Conditions (G) ETC.

(D) Type of Charger

\*5. The Centers for Medicare and Medicaid Service (CMS).  
The healthcare Common Procedure Coding System (HCPCS)

All specifications are subject to change without prior notice.

Shopperider Mobility Products, Inc. reserves the rights of any changes to the unit.

1) Include the anti-tip wheel or the rear castor.

2) Driver weight may exceed weight of the unit; speed must be reduced when turning.

3) The actual driving range varies with the factors shown as below:

- The weight of occupant
- Ground surface
- Battery capacity and conditions
- Type of charger
- Ambient temperature
- The way of driving
- If the battery and mechanical moving parts fully break in
- Etc.

# SHOPRIDER®

Toll-Free: 800.743.0772

www.shopperider.com