

SHOPRIDER®



6RUNNER® 10 DLX [858WM]

Introducing a new member to Shoprider's 6Runner® series power chairs, the 6Runner® 10 Deluxe offers the enhanced stability of all 6 points of contact with the ground and true mid-wheel drive technology. The 6Runner® 10 Deluxe beautifully captures the functionality and maneuverability for someone on the go. Ride along with confidence with the exceptional quality you can always find in a Shoprider brand power chair.

Toll-Free: 800.743.0772
www.shoprider.com



Adjustable Height



Flip-up Arm Rest



Reclining Seat

MID WHEEL DRIVE

- Luxurious leatherette seat
- Front and Rear suspension
- Stable four post seating system
- Reclining seat and adjustable seat height
- Adjustable length, height and angle foot plate
- Adjustable width, height and flip-up arm rest
- Midwheel drive maneuverability
- 10" wheels
- 300 lb. capacity

6RUNNER[®] 10 DELUXE



Red



Blue

Overall Dimensions (L x W x H)	in	36 x 22 x 38.6
Number, Size of Tires (front; rear)	Front (rear)	2, 6"; 2, 10"; 2, 5"
Suggested User Weight (on a level road)	lb	300
Total Weight (w/battery. 50 Ah*2)	lb	205
Battery Module Capacity	_V_Ah x Pcs	12V50Ah x 2
Charger	Amp	Off Board (4A)
Output Power of Motor	hp	0.8 x2
Maximum Speed ²	mph	5
Range (per charge) (After the battery and mechanical moving parts fully break in)	mile	20
Turning Radius	in	18

All specifications are subject to change without prior notice. Shoprider Mobility Products, Inc. reserves the rights of any changes to the unit.

1) Include the anti-tip wheel or the rear castor.

2) Driver weight may exceed weight of the unit; speed must be reduced when turning.

3) The actual driving range varies with the factors shown as below:

a) The weight of occupant

b) Ground surface

c) Battery capacity and conditions

d) Type of charger

e) Ambient temperature

f) The way of driving

g) If the battery and mechanical moving parts fully break in

h) Etc.