

**SHOPRIDER**



# ECHO [SL7-3]

The extremely lightweight frame of the Echo scooter allows for the convenience of taking it everywhere you want to go. The Echo is compact, affordable mobility that can take you through small spaces with ease and stability. With its quick connect battery pack, charging your scooter has never been so easy.

Toll-Free: 800.743.0772  
[www.shoprider.com](http://www.shoprider.com)



Fitted Bag in the Seat



Back cushion w/pouch, Armrest

## COMPACT SCOOTER

- 250 lb. weight capacity
- Quick connect battery pack
- Charge battery pack at any convenient location on/off the scooter
- Puncture proof tires
- Adjustable seat height
- Front Basket included

# ECHO



Overall Dimensions (L x W x H)	in	37 x 21 x 36
Number, Size of Tires (front; rear)	Front (rear)	1, 8" (2, 8")
Suggested User Weight (on a level road)	lb	250
Weight of Heaviest Part (w/o batteries, baskets and seat)	lb	35
Total Weight (w/ batteries, basket, seat)	lb	70
Each Battery Module Weight	lb	15
Battery Module Capacity	_V_Ah x Pcs	12V10Ah x 2
Charger		Off Board (2A)
Output Power of Motor	hp	0.45
Maximum Speed <sup>2</sup>	mph	3.8
Max. Climbing Angle/Safe Climbing Angle (based on user weight of 90kg/200lb.)	% (°)	15 (9) 130lbs/10 (6) 187 lbs
Range (per charge)(After the battery and mechanical moving parts fully break in)	mile	8
Turning Radius	in	32
Handle Bar	Type	T-Bar
Drive Wheel	Type	Rear Wheel

All specifications are subject to change without prior notice. Shoprider Mobility Products, Inc. reserves the rights of any changes to the unit.

1) Include the anti-tip wheel or the rear castor.

2) Driver weight may exceed weight of the unit; speed must be reduced when turning.

3) The actual driving range varies with the factors shown as below:

a) The weight of occupant

d) Type of charger

g) If the battery and mechanical moving parts fully break in

b) Ground surface

e) Ambient temperature

c) Battery capacity and conditions

f) The way of driving

h) Etc.