

Jetstream L (888WAL)

SHOPRIDER[®]



Shoprider is proud to present the new, revolutionary **Jetstream L**. With exceptional comfort, **Shoprider's Jetstream L** power chair has quiet power and maneuverability and is easily manageable with a programmable joystick controller. Outstanding stability and style make this power chair a cut above.

Available Colors



Red



Blue



Green



Titanium

On-board Charger

Footrest
Height/Length
Adjustment

Adjustable Headrest

Reclining Seat

Armrest
Height/Width
Adjustment

14" Wheel



Accessories



Cane / Oxygen Holder



Rear Basket



Walker Holder



Elevating Leg Rest (ELR)



Splash N' Dash™ Quick
Color Change

Jetstream M (888WAM) Specifications

Overall Dimensions (L' x W x H)	in	40 x 26 x 40
Number, Size of Tires (front; rear)	ea, in	2, 6; 2, 14
Weight Capacity	lb	Max. 300
Weight of Heaviest Part (without batteries)	lb	99
Total Unit Weight	lb	215
Battery Module Weight	lb	25 each x 2
Battery Module Capacity	_V_Ah x Pcs	12V35Ah x 2
Charger		On-board (4A)
Output Power of Motor	hp	1.34
Turning Radius	in	26
Ground Clearance	in	3.5
Seat to Floor Height	in	22
Maximum Speed	mph	6

All specifications are subject to change without prior notice. Shoprider Mobility Products, Inc. reserves the rights of any changes to the unit.

- 1) Include the anti-tip wheel or the rear castor.
- 2) Driver weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown as below:
 - a) The weight of occupant
 - b) Ground surface
 - c) Battery capacity and conditions
 - d) Type of charger
 - e) Ambient temperature
 - f) The way of driving
 - g) If the battery and mechanical moving parts fully break in
 - h) Etc.

Jetstream L (888WAL)

- Luxurious leatherette 20" high back recline executive seat
- Adjustable length, height and angle foot plate
- Adjustable width, height and flip-up arm rest
- On-board charger with retractable cord
- Rear wheel drive for stability
- 14" wheel
- 300 lb. capacity



Adjustable height & length footplate



On-board charger with retractable cord



Flip-up seat for battery access



Reclining seat

SHOPRIDER

Toll-Free: 800.743.0772
www.shoprider.com