

SHOPRIDER®



K0827

SHOPRIDER HD [888WHD]

The Shoprider HD is the heaviest duty of our range of power chair solutions. The 500lb weight capacity Shoprider HD will deliver solid and dependable performance expected from any powerchair product with Shopriders brand name on it. Rear drive stability, extra large standard seat and footplate, full suspension, solid tires, and extra powerful 1 horsepower motors are delivered in the HD's standard package. Shopriders trusted brand name and quality are delivered at no additional expense.

Toll-Free: 800.743.0772
www.shoprider.com



HEAVY DUTY POWERCHAIR

- Powerful twin drive motors supplying 1 horsepower each
- Programmable controller 70 amp
- Heavy-duty gearboxes
- Flat Free Tires
- Arms height and width adjustable
- Swing up arms for easy transfer
- Fully adjustable Extra Large footplate – length, height, angle
- Extra Wide seat (width 22 in and depth 21 in)

SHOPRIDER HD



Overall Dimensions (L x W x H)	in	47 x 28 x 50
Number, Size of Tires (front; rear)	Front (rear)	2, 9"; 2, 10"
Suggested User Weight (on a level road)	lb	500
Total Weight (w/75Ah*2)	lb	301
Each Battery Module Weight (75Ah)	lb	57
Battery Module Capacity	_V_Ah x Pcs	12V75Ah x 2
Charger		Off Board (8A)
Output Power of Motor	hp	1 x2
Maximum Speed ²	mph	4.375
Max. Climbing Angle/Safe Climbing Angle(based on user weight of 90kg/200lb.)	% (°)	10 (6)
Range (per charge) (After the battery and mechanical moving parts fully break in)	mile	15
Turning Radius	in	37.5

All specifications are subject to change without prior notice. Shoprider Mobility Products, Inc. reserves the rights of any changes to the unit.

1) Include the anti-tip wheel or the rear castor.

2) Driver weight may exceed weight of the unit; speed must be reduced when turning.

3) The actual driving range varies with the factors shown as below:

a) The weight of occupant

b) Ground surface

c) Battery capacity and conditions

d) Type of charger

e) Ambient temperature

f) The way of driving

g) If the battery and mechanical moving parts fully break in

h) Etc.