

**SHOPRIDER**



**K0899**

## **SPRINTER XL4 DELUXE** [889XLSBN]

The fully equipped Sprinter™ XL4 Deluxe offers many luxurious features such as a full light package, delta style tiller, and a comfy ride with its full suspension and 4 wheel design. The Sprinter™ XL4 allows for a smooth and steady ride with lots of flair for the larger rider.

Toll-Free: 800.743.0772  
[www.shoprider.com](http://www.shoprider.com)



## HEAVY DUTY SCOOTER

- 500 lb. weight capacity
- Full light package
- 12" wheels
- Up to 25 miles range
- Emergency hand brake
- 20" deluxe captain seat that slides, swivels and reclines
- 2 HP motor
- Full suspension
- Delta style adjustable tiller
- 20 speed settings
- Front basket included

# SPRINTER XL4 DELUXE

Overall Dimensions (L x W x H)	in	56 x 28 x 51
Number, Size of Tires (front; rear)	Front (rear)	2, 12" (2, 12")
Suggested User Weight (on a level road)	lb	Max. 500
Weight of Heaviest Part (w/o baskets and seat)	lb	115
Total Weight (w/ batteries, basket, seat)	lb	240
Each Battery Module Weight	lb	50
Battery Module Capacity	_V_Ah x Pcs	12V33-75Ah x 2
Charger		Off Board (8A)
Output Power of Motor	hp	2
Maximum Speed <sup>2</sup>	mph	5
Max. Climbing Angle/Safe Climbing Angle (based on user weight of 90kg/200lb.)	% (°)	14 (8)
Range (per charge) (After the battery and mechanical moving parts fully break in)	mile	25
Turning Radius	in	72
Ground Clearance	in	4
Handle Bar	Type	Delta

All specifications are subject to change without prior notice. Shoprider Mobility Products, Inc. reserves the rights of any changes to the unit.

1) Include the anti-tip wheel or the rear castor.

2) Driver weight may exceed weight of the unit; speed must be reduced when turning.

3) The actual driving range varies with the factors shown as below:

a) The weight of occupant

b) Ground surface

c) Battery capacity and conditions

d) Type of charger

e) Ambient temperature

f) The way of driving

g) If the battery and mechanical moving parts fully break in

h) Etc.