

# SUNRUNNER™ 3 DELUXE 888-3L

The Sunrunner™ 3 Deluxe is truly cutting-edge. With the vespa style tiller cover, black leatherette pillow top seat, and articulating front wheels allows you to drive with ease and luxury. The Sunrunner™ 3 Deluxe comes with the world patented on-board charger with a retractable cord. You never have to worry about leaving your charger at home. The articulating front wheel design offers four-wheel stability with three-wheel maneuverability.

## Available Colors



Burgundy



- Standard reclining and sliding seat
- 300 lb. weight capacity
- Full light package
- On-board retractable charger
- Adjustable length
- Full light package including blinkers/hazards
- Seat height adjustable
- C-bar style handlebar
- Adjustable tiller
- Basket included

Overall Dimensions (L1 x W x H)	in	48-51 x 24 x 41
Number, Size of Tires (front, rear)	Front (Rear)	2, 9" (2, 10")
Suggested User Weight (on a level road)	lb (kg)	Max. 300 (136)
Weight of Heaviest Part (w/o baskets and seat)	lb (kg)	N/A
Total Weight (w/ batteries, basket, seat)	lb (kg)	160 (73)
Each Battery Module Weight	lb (kg)	24 (11)
Battery Module Capacity	_V_Ah x Pcs	12V36Ah x 2
Charger		On Board (4A)
Output Power of Motor	hp (W)	.80 (597)
Maximum Speed <sup>2</sup>	mph (km/hr)	5 (8)
Maximum Climbing Angle / Safe Climbing Angle(based on user weight of 90kg/200lb)	% (°)	14 (8) / 21 (12)*
Range (per charge)(After the battery and mechanical moving parts fully break in)	mile (km)	24 (38)
Turning Radius	in (mm)	39 (991)
Ground Clearance	in (mm)	3 (76)
Handle Bar	Type	C-Bar

**SHOPRIDER**

Toll-Free: 800.743.0772

www.shoprider.com