

**SHOPRIDER®**



**K0823**

# STREAMER SPORT [888WA]

The Streamer Sport comes with puncture proof 10" tires, on board charger and a dramatically sloped shroud. With only the best technology, the Streamer Sport is equipped with self diagnosing P&G VSI controllers. The Streamer Sport is full of features.

Toll-Free: 800.743.0772  
[www.shoprider.com](http://www.shoprider.com)



Pan Seat Option

## REAR WHEEL DRIVE

- Adjustable angle foot plate
- Midwheel maneuverability
- 10" wheel
- 300 lb. capacity
- Height and angle adjustment arm rests
- Optional Swing away controller arm bracket

# STREAMER SPORT



Red

Overall Dimensions (L x W x H)	in	36 x 24 x 45
Number, Size of Tires (front; rear)	Front (rear)	2, 6"; 2, 10"
Suggested User Weight (on a level road)	lb	300
Total Weight (w/36 Ah batteries, seat)	lb	203
Each Battery Module Weight	lb	5.8
Battery Module Capacity	_V_Ah x Pcs	12V35Ah x 2
Charger		Off Board (5A)
Output Power of Motor	hp	0.60 x2
Maximum Speed <sup>2</sup>	mph	5
Max. Climbing Angle/Safe Climbing Angle(based on user weight of 90kg/200lb.)	% (°)	10 (6)
Range (per charge) (After the battery and mechanical moving parts fully break in)	mile	20
Turning Radius	in	17.7

All specifications are subject to change without prior notice. Shoprider Mobility Products, Inc. reserves the rights of any changes to the unit.

1) Include the anti-tip wheel or the rear castor.

2) Driver weight may exceed weight of the unit; speed must be reduced when turning.

3) The actual driving range varies with the factors shown as below:

a) The weight of occupant

b) Ground surface

c) Battery capacity and conditions

d) Type of charger

e) Ambient temperature

f) The way of driving

g) If the battery and mechanical moving parts fully break in

h) Etc.