

SHOPRIDER®




K0822

STREAMER SPORT REHAB [888WA-R]

The Streamer Sport Rehab comes with puncture proof 10" tires, and a fully adjustable rehab seating with pan seat. With only the best technology, the Streamer Sport Rehab is equipped with self diagnosing PG Drives VSI controllers. The Streamer Sport Rehab is full of features.

Toll-Free: 800.743.0772
www.shoprider.com




Seat Cushion Option

REAR WHEEL DRIVE

- 300 lb. capacity
- 10" puncture proof tires
- Angle adjustment footplate
- Height, length and angle adjustment arm rests
- Swing away controller arm bracket
- Height adjustable seat
- Fully adjustable Rehab seat

STREAMER SPORT REHAB



Red

Blue

Overall Dimensions (L x W x H)	in	37 x 24 x 32
Number, Size of Tires (front; rear)	Front (rear)	2, 6"; 2, 10"
Suggested User Weight (on a level road)	lb	300
Total Weight (w/36 Ah batteries, seat)	lb	180
Each Battery Module Weight	lb	4.5-5
Battery Module Capacity	_V_Ah x Pcs	12V35Ah x 2
Charger		On Board (3A) or Off board (5A)
Output Power of Motor	hp	0.60 x2
Maximum Speed ²	mph	5
Max. Climbing Angle/Safe Climbing Angle(based on user weight of 90kg/200lb.)	% (°)	10 (6)
Range (per charge) (After the battery and mechanical moving parts fully break in)	mile	20
Turning Radius	in	17.7

All specifications are subject to change without prior notice. Shoprider Mobility Products, Inc. reserves the rights of any changes to the unit.

1) Include the anti-tip wheel or the rear castor.

2) Driver weight may exceed weight of the unit; speed must be reduced when turning.

3) The actual driving range varies with the factors shown as below:

a) The weight of occupant

b) Ground surface

c) Battery capacity and conditions

d) Type of charger

e) Ambient temperature

f) The way of driving

g) If the battery and mechanical moving parts fully break in

h) Etc.