

STREAMER SPORT REHAB (888WA-R) K0822

The Streamer Sport Rehab comes with puncture proof 10" tires, and a fully adjustable rehab seating with pan seat. With only the best technology, the Streamer Sport Rehab is equipped with self diagnosing PG Drives VSI controllers. The Streamer Sport Rehab is full of features.

Available Colors



Red



Blue



Seat Cushion Option shown

- 300 lbs. weight capacity
- 10" puncture proof tires
- Available in Red and blue
- Angle adjustment footplate
- Height, length and angle adjustment arm rests
- Swing away controller arm bracket
- Height adjustable seat
- Fully adjustable Rehab seat

Overall Dimensions (L ¹ x W x H ²)	in	35 x 24 x 35.4
Number, Size of Tires (front; rear)	Front (Caster)/Middle (Driver)	2, 6" (2, 10")
Total Weight (w/36Ah batteries, basket, seat)	lbs	180 lbs
Battery Capacity	_V_Ah x Pcs	12V35Ah x2
Charger	Amp	3
Output Power of Motor (Reference only)	hp	0.6 x2
Max. Single User Weight (on the level road)	lb	300
Maximum Speed ³ (on the level road)	mph	4.5 - 5
Turning Radius	in	17.7
Range ⁴ (per charge with 36Ah Battery / 90 Kg User weight) (After the battery and mechanical moving parts fully break in)	mile	up to 20 miles

All specification subject to change without prior notice. Shoprider Mobility Products, Inc. reserves the rights of any changes to the unit.

1) Include footplate & anti-tip wheel.

2) Depending on the height of the seat back.

3) When proceeding up any incline, please move your seat to the most forward position or if you have the deluxe seat, make sure that it is in the 90 degree (upright) position.

4) The actual driving range varies with the factors shown as below:

a) The weight of occupant
b) Level ground
c) Battery conditions

d) Type of charger
e) Ambient temperature
f) The way of driving

g) Etc.

SHOPRIDER

Toll-Free: 800.743.0772

www.shoprider.com