

**Wizz** (888WNL)

**SHOPRIDER**



The new, redesigned **Wizz** Mid-wheel Power Chair offers great comfort and style while providing unmatched maneuverability. Its true Mid-wheel design and **Shoprider** value provide the best overall powerchair package in today's industry. The larger and luxuriously comfortable pillow-top reclining seat and various accessory options vault the **Wizz** into a category of its own.

## Available Colors



Red



Blue



**The Wizz fits  
your lifestyle!**



Adjustable Headrest

Reclining  
Seat

Armrest Height  
and Length  
Adjustment

On-board Charger

Adjustable  
Length Footrest

10" Wheel

## Accessories



Cane / Oxygen Holder



Rear Basket



Elevating Leg Rest (ELR)



Walker Holder



Adjustable length  
footplate



On-board charger



Swing away controller



Reclining seat

## Wizz (888WNL) Specifications

Overall Dimensions (L' x W x H)	in	41 x 23 x 40
Number, Size of Tires (front; mid; rear)	ea, in	2, 6; 2, 10; 2, 5
Weight Capacity	lb	Max. 300
Weight of Heaviest Part	lb	99
Total Unit Weight	lb	203
Battery Module Weight	lb	25 each x 2
Battery Module Capacity	_V_Ah x Pcs	12V35Ah x 2
Charger		On Board (3A)
Output Power of Motor	hp	.75
Turning Radius	in	21
Ground Clearance	in	3.5
Seat to Floor Height	in	21
Maximum Speed	mph	5

All specifications are subject to change without prior notice. Shoprider Mobility Products, Inc. reserves the rights of any changes to the unit.

1) Include the anti-tip wheel or the rear castor.  
2) Driver weight may exceed weight of the unit; speed must be reduced when turning.  
3) The actual driving range varies with the factors shown as below:  
a) The weight of occupant  
b) Ground surface  
c) Battery capacity and conditions  
d) Type of charger  
e) Ambient temperature  
f) The way of driving  
g) If the battery and mechanical moving parts fully break in  
h) Etc.

## Wizz (888WNL)

- Luxurious leatherette 20" high back recline executive seat
- Adjustable length and angle foot plate
- Adjustable width and length arm rest
- Mid-wheel drive
- 10" mid-wheel
- 300 lb. capacity

**SHOPRIDER**

Toll-Free: 800.743.0772

www.shoprider.com