

# PRICE LIST & ORDER FORM



For easy ordering and information, please contact Customer Service Toll Free at 1.800.743.0772 FAX: 1.877.797.7081 Effective July 2009



**UL8W SMARTIE**

NAME: \_\_\_\_\_  
 PHONE NO: \_\_\_\_\_  
 ACCT NO: \_\_\_\_\_  
 PO #: \_\_\_\_\_ DATE: \_\_\_\_\_  
 ADDRESS: \_\_\_\_\_  
 CITY: \_\_\_\_\_ STATE: \_\_\_\_\_ ZIP: \_\_\_\_\_  
 SHIP TO: \_\_\_\_\_  
 PHONE NO: \_\_\_\_\_  
 ADDRESS: \_\_\_\_\_  
 CITY: \_\_\_\_\_ STATE: \_\_\_\_\_ ZIP: \_\_\_\_\_

## SPECIFICATIONS

### OVERALL DIMENSIONS

Base Length.....36"  
 Base Width.....23"  
 Height.....37"  
 (Ground to Seat Top)  
 Ground Clearance.....2"

### WEIGHTS

User Weight Capacity.....250Lb  
 Front Chassis w/ batteries.....27.1 Lbs  
 Rear Chassis w/ drive train.....29 Lbs  
 Total Weight(w/ Bat's & Seat) .....93 Lbs  
 Battery pack weight.....5.9 Lbs  
 Captain seat w/controller assy.....37.1 Lbs  
 Battery Module Weight (Pack).....5.9 Lbs

### TIRES

Number, Size of Front Tire.....2,6"  
 Number, Size of Rear.....2,8"

### BATTERIES

Envirofriendly Battery Pack.....5.9 Lbs

### CHARGER

Charger-Off Board.....2 amp

### MOTOR POWER OUTPUT

Horsepower.....0.45 x2

### DRIVING SPECS

Turning Radius .....15.5"  
 Range .....10 mph  
 Maximum Speed .....3.75 mph  
 Maximum Climbing Angle .....10% / 6"

All specifications are subject to change without notice. Shoprider Mobility Products, Inc® reserves the rights of any changes to the unit:

- 1) Include the anti-tip wheel or rear castor.
- 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown below:
  - a) The weight of user.
  - b) Ground surface and terrain.
  - c) Battery capacity and conditions.
  - d) Type of charger.
  - e) Ambient Temperature.
  - f) The way and/or method of driving.
  - g) If the battery and mechanical moving parts are fully broken in.

### BASE SCOOTER

UL8W

Smartie Powerchair  
 Red or Blue

K0821

### MSRP

\$5,495.00

### BATTERY

PH012010D01

Envirofriendly Battery Pk

x1

\$740.00\*\*\*

\*\*\*Envirofriendly Battery must be purchased with Smartie

### SEATS (W"xD"xH")

SV1817½19G

18" Grey Vinyl Deluxe Captain Seat

\$400.00\*\*

### ARM REST ASSY

104300-00106-A3

(L) Arm Rest Assy w/ Arm Pad

\$140.00\*\*

104400-00106-A3

(R) Arm Rest Assy w/ Arm Pad

\$140.00\*\*

### TIRES

105201-69101

200X70 Solid Rear Wheel Assy

x2

\$90.00\*\*ea

1351-00900-00

6" Solid Tire Front

x2

\$140.00\*\*ea

### ANTI-TIP WHEELS

209928-89300

Anti-Tip Wheel Assy

x2

\$20.00\*\*ea

### CHARGER

109201-88805

3 pin; 2 Amp charger; UL8W; Li

\$100.00\*\*

### CONTROLLER

1361-00300-01

Joystick Controller-UL8WPB

\$998.00\*\*

### FOOT REST

2329-00100-11

Full Function Foot Plate

\$200.00\*\*

ELRUL-R(K)

Elevating Leg Rest Kit in Red

\$335.00

ELRUL-B(K)

Elevating Leg Rest Kit in Black

\$335.00

### ACCESSORIES

OXY-03

Mesh Oxygen Tank Holder

\$60.00

CCH-003-GK(K)

Cane Holder Kit

\$150.00

WH-0185-GK(K)

Walker Holder Kit

\$150.00

\*\*STANDARD WITH UNIT

TOTAL \$ \_\_\_\_\_