



PRICE LIST & ORDER FORM

Effective November 2006

For easy ordering and information contact Customer Service Toll Free at:

800-743-0772 Fax: 877-797-7081

| | |
|----------------|------------------|
| Name _____ | Phone No _____ |
| Acct No _____ | PO # _____ |
| Address _____ | City/State _____ |
| Ship To _____ | Phone No _____ |
| Address _____ | City/State _____ |
| Comments _____ | PO # _____ |

Stock For your convenience, order by stock number (e.g. BAT-U1)

BASE CHAIR

- 888WAM11-BGRD JETSTREAM-M K0849 \$5,395.00
- 888WAM11-BLUE JETSTREAM-M K0849 \$5,395.00

BATTERY

- BAT-U1 with seat pan upgrade K0848

- 12V 35AMP x2 \$220.00

SEATS (W"xD"xH")

- SV201820RB 20x18x20 Recline High-Back Pillow-Top \$370.00**
 - SC201822RG 20x18 Recline High-Back Grey Cloth \$150.00
 - SP202023RG 20w x 20d Seat Pan/Capt back-Recline \$150.00
- also available in 16" & 18" *Seat Cushion not included

ARM REST ASSY

- 104400-00100-B2 (R) Flipback Adj Arm Rest Assy \$160.00**
- 104300-00100-B2 (L) Flipback Adj Arm Rest Assy \$160.00**

TIRES

- TSR26085 10" Flat-Free Rear Tires x2 \$100.00**

CASTORS

- 1351-00601-00 6" Front Castor Assy x2 \$280.00**

ANTI-TIP WHEELS

- 800405-002-01 Anti-Tip Wheel Assy x2 Please Call

FOOT REST

- Please Call Flip Up Foot Plate Standard
- ELRJET-B(K) Elevating Leg Rest Kit in Black \$232.00
- ELRJET-R(K) Elevating Leg Rest Kit in Red \$232.00
- 202934-00104 + (L) Leg Rest Hanger-ELR/SAFR \$32.00
- 5008-08020-210 Hanger Tightening Knob x2 \$24.00
- 202934-00105 + (R) Leg Rest Hanger-ELR/SAFR \$32.00
- 5008-08020-210 Hanger Tightening Knob x2 \$24.00
- 800907-007-00 Swing Away Foot Rest (R) \$80.00
- 800907-007-01 Swing Away Foot Rest (L) \$80.00

CONTROLLER

- 1361-00111-00A2 VSI Joystick Controller \$998.00**
- 1343-PHFW-SQ VSI Controller Bracket \$36.00**

CHARGER

- 109201-88115 4 amp Charger On-Board \$200.00**

ACCESSORIES

- CCH-003(K) Cane Holder \$150.00
- CH-02(K) Oxygen Holder (Metal) \$150.00
- WH-0185(K) Walker Holder \$150.00
- OXY-03(K) Mesh Oxygen Tank Holder \$60.00
- 300909-03 Safety Belt \$34.00**
- 19903-00302 Front Castor Curb Guard-Set \$70.00

TOTAL \$

**STANDARD WITH UNIT

JETSTREAM-M SPECIFICATIONS

OVERALL DIMENSIONS

- Base Length.....40"
- Base Width.....26"
- Height.....36"
- (Ground to Seat Top)
- Ground Clearance.....2"

WEIGHTS

- User Weight Capacity.....300 Lbs
- Total Unit Weight.....205 Lbs
- Weight of Heaviest Part.....89 Lbs
- Battery Module Weight (U1).....25 Lbs

TIRES/CASTORS

- Front Castor Assy.....6"
- Rear Flat-Free Tires.....10"

BATTERIES

- BAT-U1 Module Capacity.....12v35Ah x2

CHARGER

- Charger-On Board.....4 amp

MOTOR POWER OUTPUT

- Horsepower.....0.8 x2

DRIVING SPECS

- Turning Radius.....24"
- Range.....up to 19 miles
- Maximum Speed.....5 mph
- Maximum Climbing Angle.....10% / 6°

All specifications are subject to change without notice. Shoprider Mobility Products, Inc® reserves the rights of any changes to the unit:

- 1) Include the anti-tip wheel or rear castor.
- 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown below:
 - a) The weight of user.
 - b) Ground surface and terrain.
 - c) Battery capacity and conditions.
 - d) Type of charger.
 - e) Ambient Temperature.
 - f) The way and/or method of driving.
 - g) If the battery and mechanical moving parts are fully broken in.