



PRICE LIST & ORDER FORM Effective November 2006

For easy ordering and information contact Customer Service Toll Free at:
800-743-0772 Fax: 877-797-7081

Name _____	Phone No _____
Acct No _____ PO # _____	Date _____
Address _____	City/State _____ Zip _____
Ship To _____	Phone No _____
Address _____	City/State _____ Zip _____
Comments _____	PO # _____

Stock For your convenience, order by stock number (e.g. BAT-75)

BASE CHAIR

<input type="checkbox"/>	888WNLLHD	6 RUNNER 14" Red	K0825	\$6,595.00
		with seat pan upgrade	K0824	add \$75.00

BATTERY

<input type="checkbox"/>	BAT-75	12V 75Ah	x2	\$250.00 ea
<input type="checkbox"/>	BAT-50	12V 50Ah	x2	\$200.00 ea

SEATS (W"xD"xH")

<input type="checkbox"/>	SV212121½RG	21Wx21D Grey Vinyl Recline Capt Seat		\$400.00 **
<input type="checkbox"/>	SV191820RG	Upgrade 19x18 Recline High-Back Grey Vinyl		\$150.00
<input type="checkbox"/>	SC201822RG	Upgrade 20x18 Recline High-Back Grey Cloth		\$150.00
<input type="checkbox"/>	SP181821RG	Upgrade 18"Seat Pan/Capt back/rcing*		\$150.00

*Also available in 16" & 20" cushion not included

ARM REST ASSY

<input type="checkbox"/>	104400-00100-B2	(R) Flipback Adj Arm Rest Assy		\$160.00 **
<input type="checkbox"/>	104300-00100-B2	(L) Flipback Adj Arm Rest Assy		\$160.00 **

TIRES

<input type="checkbox"/>	TSFR3008G Grey	14" Flat-Free Tires Mid(3.00x8)x2		\$112.00 ea **
--------------------------	----------------	-----------------------------------	--	----------------

CASTORS

<input type="checkbox"/>	105101-69000	7" Front Castor Assy	x2	\$120.00 ea **
<input type="checkbox"/>	105101-69000	7" Rear Castor Assy	x2	\$120.00 ea **

FOOT REST

<input type="checkbox"/>	2329-00100-10	Flip Up Foot Plate		\$200.00 **
<input type="checkbox"/>	ELRJET-R(K)	Elevating Leg Rest Kit in Red		\$232.00
<input type="checkbox"/>	ELRJET-B(K)	Elevating Leg Rest Kit in Black		\$232.00
<input type="checkbox"/>	202934-00104 +	(L) Leg Rest Hanger-ELR/SAFR		\$32.00
<input type="checkbox"/>	202934-00105 +	(R) Leg Rest Hanger-ELR/SAFR		\$32.00
<input type="checkbox"/>	800907-007-00	Swing Away Foot Rest (R)red only		\$80.00
<input type="checkbox"/>	800907-007-01	Swing Away Foot Rest (L)red only		\$80.00

CONTROLLER

<input type="checkbox"/>	1361-00111-00A2	VSI Joystick Controller		\$998.00 **
<input type="checkbox"/>	1343-PHFW-SQ	VSI Controller Bracket		\$36.00 **

CHARGER

<input type="checkbox"/>	109201-88500	8 amp Charger Off-Board		\$300.00 **
--------------------------	--------------	-------------------------	--	-------------

ACCESSORIES

<input type="checkbox"/>	CCH-003(K)	Cane Holder		\$150.00
<input type="checkbox"/>	CH-02(K)	Oxygen Tank Holder (Metal)		\$150.00
<input type="checkbox"/>	WH-0185(K)	Walker Holder		\$150.00
<input type="checkbox"/>	OXY-03(K)	Mesh Oxygen Tank Holder		\$60.00
<input type="checkbox"/>	300909-03	Safety Belt		\$34.00 **
<input type="checkbox"/>	RB-0125(K)	Rear Basket		\$170.00

**STANDARD WITH UNIT

TOTAL \$

6 RUNNER SPECIFICATIONS

OVERALL DIMENSIONS 888WNLLHD

Base Length.....	41"
Base Width.....	25"
Height.....	21", 22", 23"
(Ground to Top of Seat)	
Ground Clearance.....	2.75"

WEIGHTS

User Weight Capacity.....	450 Lbs
Total Unit Weight w/ 75Ah	277 Lbs
Weight of Heaviest Part.....	82 Lbs
Battery Module Weight (75Ah)....	57 Lbs ea

TIRES/CASTORS

Front Castor Assy.....	7"
Mid Wheel Tire Flat Free.....	14"
Rear Castor Assy.....	7"

BATTERIES

75Ah Module Capacity.....	
50Ah Module Capacity.....	

CHARGER

Charger Off Board.....	8 amp
------------------------	-------

MOTOR POWER OUTPUT

Horsepower.....	1.34 x2
-----------------	---------

DRIVING SPECS

Turning Radius.....	25"
Range w/ 75Ah batteries.....	up to 20 mi
Maximum Speed.....	3.75 mph
Maximum Climbing Angle.....	10% / 6°

All specifications are subject to change without notice. Shoprider Mobility Products, Inc® reserves the rights of any changes to the unit:

1) Include the anti-tip wheel or rear castor.
 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.

3) The actual driving range varies with the factors shown below:

- a) The weight of user.
- b) Ground surface and terrain.
- c) Battery capacity and conditions.
- d) Type of charger.
- e) Ambient Temperature.
- f) The way and/or method of driving.
- g) If the battery and mechanical moving parts are fully broken in.