

PRICE LIST & ORDER FORM



For easy ordering and information, please contact Customer Service **Toll Free at 1.800.743.0772 FAX: 1.877.797.7081** Effective July 2009



WNLLHD 6RUNNER14

NAME: _____
 PHONE NO: _____
 ACCT NO: _____
 PO #: _____ DATE: _____
 ADDRESS: _____
 CITY: _____ STATE: _____ ZIP: _____
 SHIP TO: _____
 PHONE NO: _____
 ADDRESS: _____
 CITY: _____ STATE: _____ ZIP: _____

SPECIFICATIONS

OVERALL DIMENSIONS

Base Length.....42"
 Base Width.....26"
 Height.....49"
 Ground to Seat Top
 Ground Clearance.....2.75"

WEIGHTS

User Weight Capacity.....450 Lbs
 Total Unit Weight(w/ 75Ah).....277 Lbs
 Weight of Heaviest Part.....82 Lbs
 Battery Weight Module (75Ah).....57 Lbs ea

TIRES/ CASTORS

Front castor Assy.....7"
 Mid Wheel Tire Flat Free.....14"
 Rear Castor Assy.....7"

BATTERIES

50Ah Module Capacity.....12v50Ah x2
 75Ah Module Capacity.....12v75Ah x2

CHARGER

Charger-Off Board.....5 amp

MOTOR POWER OUTPUT

Horsepower.....1.34 x2

DRIVING SPECS

Turning Radius23.5"
 Range.....up to 26 mph
 Maximum Speed.....3.75 mph
 Maximum Climbing Angle.....10% / 6°

BASE SCOOTER

888WNLLHD 6 RUNNER 14" Red K0825 \$6,595.00
 with seat pan upgrade K0824 \$6,745.00

BATTERY

BAT-50 12V 50Ah x2 \$200.00 ea**
 BAT-75 12V 75Ah x2 \$300.00 ea

SEATS (W"xD"xH")

SV212121½RG 21Wx21D Grey Vinyl Recline Capt Seat \$400.00 **
 SP161618RG Upgrade 16" Seat Pan/Capt back/rcing* \$400.00
 SP181821RG Upgrade 18" Seat Pan/Capt back/rcing* \$400.00
 SP202023RG Upgrade 20" Seat Pan/Capt back/rcing* \$400.00
 SV201820RB 20x18x20 Recline Black Pillow Top \$370.00
 SV221922RG 22x19x22 Grey Vinyl Recline Capt Seat \$390.00
 SV191820RG Upgrade 19x18 Recline High-Back Grey Vinyl \$400.00
 SV212121½RB 21Wx21D Black Vinyl Recline Capt Seat \$400.00
 SV1817½19G 18x17½x19 Recline Capt. Seat- Grey Vinyl \$400.00
 SV2319½19½B 23x19½x19½ Recline Capt. Seat- Black Vinyl \$420.00
 *Cushions not included

ARM REST ASSY

104300-00112 (L) Flipback Adj Arm Rest Assy \$140.00**
 104400-00112 (R) Flipback Adj Arm Rest Assy \$140.00**

TIRES

TSFR3008G Grey 14" Flat-Free Tires Mid(3.00x8) x2 \$112.00ea **

CASTORS

105101-69000 7" Front Castor Assy x2 \$60.00 ea **
 105101-69000 7" Rear Castor Assy x2 \$60.00 ea **

FOOT REST

2329-00100-10 Flip Up Foot Plate \$200.00 **
 ELRJET-R(K) Elevating Leg Rest Kit in Red \$232.00
 ELRJET-B(K) Elevating Leg Rest Kit in Black \$232.00
 202934-00104 + (L) Leg Rest Hanger-ELR/SAFR \$32.00
 202934-00105 + (R) Leg Rest Hanger-ELR/SAFR \$32.00
 800907-007-00 Swing Away Foot Rest (R)red only \$80.00
 800907-007-01 Swing Away Foot Rest (L)red only \$80.00

CONTROLLER

1361-00114-00 VR2 Controller \$360.00**
 1361-00204-02 VR2 Control Module 60amp \$776.00**

CHARGER

109201-88124 5 amp Charger Off-Board \$230.00 **

ACCESSORIES

CCH-003(K) Cane Holder \$150.00
 0216(K) Oxygen Tank Holder (Metal) \$150.00
 WH-0185(K) Walker Holder \$150.00
 OXY-03(K) Mesh Oxygen Tank Holder \$60.00
 300909-03 Safety Belt \$34.00**
 RB-0125(k) Rear Basket 20"L x 14½"W x 8"D \$170.00

All specifications are subject to change without notice. Shoprider Mobility Products, Inc.® reserves the rights of any changes to the unit:

- 1) Include the anti-tip wheel or rear castor.
- 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown below:
 - a) The weight of user.
 - b) Ground surface and terrain.
 - c) Battery capacity and conditions.
 - d) Type of charger.
 - e) Ambient Temperature.
 - f) The way and/or method of driving.
 - g) If the battery and mechanical moving parts are fully broken in.

**STANDARD WITH UNIT

TOTAL \$ _____