



For easy ordering and information contact Customer Service Toll Free at:
800-743-0772 Fax: 877-797-7081

Name _____	Phone No _____
Acct No _____	PO # _____
Date _____	
Address _____	City/State _____ Zip _____
Ship To _____	Phone No _____
Address _____	City/State _____ Zip _____
Comments _____	PO # _____

Stock For your convenience, order by stock number (e.g. BAT-U1)

BASE CHAIR

- 888WNLM-BGRD 6 RUNNER 10" Burgundy Red K0823 \$5,595.00
- 888WNLM-BLUE 6 RUNNER 10" Blue K0823 \$5,595.00

BATTERY

- BAT-U1 with seat pan upgrade K0822 x2 \$220.00

SEATS (W"xD"xH")

- SV201820RB 20X18X20 Recline Blk Pillow Top Standard
- SV2121211/2RG 21x21x211/2 Recline Capt Seat-Grey vinyl \$200.00
- SV2121211/2RB 21x21x211/2 Recline Capt Seat-Black vinyl \$200.00
- SV191820RG Upgrade 19x18 Recline High-Back Grey Vinyl \$150.00
- SC201822RG Upgrade 20x18 Recline High-Back Grey Cloth \$150.00
- SP181821RG Upgrade 18"Seat Pan/Capt back/reclining \$150.00

ARM REST ASSY

- 104400-00100-B2 (L) Flipback Adjustable Arm Rest Assy \$160.00**
- 104300-00100-B2 (R) Flipback Adjustable Arm Rest Assy \$160.00**

TIRES

- TSR26085 Grey 10" Flat-Free Drive Tires(260x85) x2 \$100.00**

CASTORS

- 1351-00800-01 7" Front Castor Assy-solid x2 \$280.00**

REAR CASTORS

- 1351-00800-01 7" Rear Castor Assy-solid x2 \$280.00**

FOOT REST

- 2329-00100-10 Flip Up Foot Plate Standard
- ELRJET-R(K) Elevating Leg Rest Kit-Red \$232.00
- ELRJET-B(K) Elevating Leg Rest Kit-Black \$232.00
- 202934-00104 + (L) Leg Rest Hanger-ELR/SAFR \$32.00
- 202934-00105 + (R) Leg Rest Hanger-ELR/SAFR \$32.00
- 800907-007-00 Swing Away Foot Rest (R) \$80.00
- 800907-007-01 Swing Away Foot Rest (L) \$80.00

CONTROLLER

- 1361-00111-00A2 Integral VSI Joystick Controller \$998.00**
- 1343-PHFW-SQ VSI Controller Bracket \$36.00**

CHARGER

- 109201-88109-A2 4 amp Charger Off-Board \$200.00**

ACCESSORIES

- CCH-003(K) Cane Holder \$150.00
- CH-02(K) Oxygen Tank Holder (Metal) \$150.00
- WH-0185(K) Walker Holder \$150.00
- OXY-03(K) Mesh Oxygen Tank Holder \$60.00
- 300909-03 Safety Belt \$34.00**
- RB-0125(k) Rear Basket 20"L x 14½"W x 8"D \$170.00

TOTAL \$

**STANDARD WITH UNIT

6 RUNNER 10 SPECIFICATIONS

OVERALL DIMENSIONS

Base Length.....30"
 Base Width.....23"
 Height.....44"
 (Ground to Seat Top)
 Ground Clearance.....3"

STANDARD SEAT

Black Vinyl, semi-recline, high back, adjustable ht head rest
 Seat W x D x H.....20"x18"x20"
 Seat Pad Cushion Height.....4"

WEIGHTS

User Weight Capacity.....300 Lbs
 Total Unit Weight.....191 Lbs
 Weight w/o Frt Rigging & Batts...133 Lbs
 Battery Module Weight (U1).....25 Lbs

TIRES/CASTORS

Front Castor Assy.....7"
 Mid Wheel Tire Flat Free.....10"
 Rear Flat Free.....7"

BATTERIES

BAT-U1 Module Capacity.....12v35Ah x2

CHARGER

Charger Off Board.....4 amp

MOTOR POWER OUTPUT

Horsepower.....0.8 x2**

DRIVING SPECS

Turning Radius.....15.75"
 Range.....up to 18.8 miles
 Maximum Speed.....5 mph
 Maximum Climbing Angle.....10% / 6°

All specifications are subject to change without notice.
 Shoprider Mobility Products, Inc® reserves the rights of any changes to the unit:
 1) Include the anti-tip wheel or rear castor.

2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.

3) The actual driving range varies with the factors shown below:

- a) The weight of user.
- b) Ground surface and terrain.
- c) Battery capacity and conditions.
- d) Type of charger.
- e) Ambient Temperature.
- f) The way and/or method of driving.
- g) If the battery and mechanical moving parts are fully broken in.