



**6 RUNNER 10**

**888WNLM**

**PRICE LIST & ORDER FORM**

Effective July 2009

For easy ordering and information contact Customer Service Toll Free at:

**800-743-0772 Fax: 877-797-7081**

Name _____	Phone No _____
Acct No _____	PO # _____
Address _____	Date _____
City/State _____	Zip _____
Ship To _____	Phone No _____
Address _____	City/State _____
City/State _____	Zip _____
Comments _____	PO # _____

**Stock** For your convenience, order by stock number (e.g. BAT-U1)

**BASE CHAIR**

<input type="checkbox"/>	888WNLM-BGRD	6 RUNNER 10" Burgundy Red	K0823	\$5,595.00
<input type="checkbox"/>	888WNLM-BLUE	6 RUNNER 10" Blue	K0823	\$5,595.00

**BATTERY**

<input type="checkbox"/>	BAT-U1	with seat pan upgrade	K0822	x2	\$140.00 ea.
--------------------------	--------	-----------------------	-------	----	--------------

**SEATS (W"xD"xH")**

<input type="checkbox"/>	SV201820RB	20X18X20 Recline Blk Pillow Top		\$370.00**
<input type="checkbox"/>	SV2121211/2RG	21x21x211/2 Recline Capt Seat-Grey vinyl		\$400.00
<input type="checkbox"/>	SV2121211/2RB	21x21x211/2 Recline Capt Seat-Black vinyl		\$400.00
<input type="checkbox"/>	SV191820RG	Upgrade 19x18 Recline High-Back Grey Vinyl		\$400.00
<input type="checkbox"/>	SC201822RG	Upgrade 20x18 Recline High-Back Grey Cloth		\$390.00
<input type="checkbox"/>	SP181821RG	Upgrade 18"Seat Pan/Capt back/reclining		\$400.00

**ARM REST ASSY**

<input type="checkbox"/>	104400-00100-B2	(R) Flipback Adjustable Arm Rest Assy		\$160.00**
<input type="checkbox"/>	104300-00100-B2	(L) Flipback Adjustable Arm Rest Assy		\$160.00**
<input type="checkbox"/>	104310-00102	Arm pad w/ height-adjustment tube		\$55.00

**TIRES**

<input type="checkbox"/>	TSR26085	Grey	10" Flat-Free Drive Tires(260x85)	x2	\$50.00**ea.
--------------------------	----------	------	-----------------------------------	----	--------------

**FRONT CASTORS**

<input type="checkbox"/>	1351-00800-01	7" Front Castor Assy-solid w/fork	x2	\$140.00**ea.
--------------------------	---------------	-----------------------------------	----	---------------

**REAR CASTORS**

<input type="checkbox"/>	1351-00800-01	7" Rear Castor Assy-solid w/fork	x2	\$140.00**ea.
--------------------------	---------------	----------------------------------	----	---------------

**CASTORS ONLY**

<input type="checkbox"/>	105101-69000	7" castor only		\$60.00 ea.
--------------------------	--------------	----------------	--	-------------

**FOOT REST**

<input type="checkbox"/>	2329-00100-10	Flip Up Foot Plate		\$200.00**
<input type="checkbox"/>	ELRJET-R(K)	Elevating Leg Rest Kit-Red		\$232.00
<input type="checkbox"/>	ELRJET-B(K)	Elevating Leg Rest Kit-Black		\$232.00
<input type="checkbox"/>	202934-00104	(L) Leg Rest Hanger-ELR/SAFR		\$32.00
<input type="checkbox"/>	202934-00105	(R) Leg Rest Hanger-ELR/SAFR		\$32.00
<input type="checkbox"/>	800907-007-00	Swing Away Foot Rest (R)		\$80.00
<input type="checkbox"/>	800907-007-01	Swing Away Foot Rest (L)		\$80.00

**CONTROLLER**

<input type="checkbox"/>	1361-00111-00A2	Integral VSI Joystick Controller		\$998.00**
<input type="checkbox"/>	1343-PHFW-SQ	VSI Controller Bracket		\$36.00**

**CHARGER**

<input type="checkbox"/>	109201-88109-A2	4 amp Charger Off-Board		\$230.00**
--------------------------	-----------------	-------------------------	--	------------

**ACCESSORIES**

<input type="checkbox"/>	CCH-003(K)	Cane Holder		\$150.00
<input type="checkbox"/>	CH-02(K)	Oxygen Tank Holder (Metal)		\$150.00
<input type="checkbox"/>	WH-0185(K)	Walker Holder		\$150.00
<input type="checkbox"/>	OXY-03(K)	Mesh Oxygen Tank Holder		\$60.00
<input type="checkbox"/>	300909-03	Safety Belt		\$34.00**
<input type="checkbox"/>	RB-0125(k)	Rear Basket 20"L x 14½"W x 8"D		\$170.00

**6 RUNNER 10 SPECIFICATIONS**

**OVERALL DIMENSIONS**

Base Length.....	30"
Base Width.....	23"
Height.....	44"
(Ground to Seat Top)	
Ground Clearance.....	3"

**STANDARD SEAT**

Black Vinyl, semi-recline, high back, adjustable ht head rest	
Seat W x D x H.....	20"x18"x20"
Seat Pad Cushion Height.....	4"

**WEIGHTS**

User Weight Capacity.....	300 Lbs
Total Unit Weight.....	191 Lbs
Weight w/o Frt Rigging & Batts...	133 Lbs
Battery Module Weight (U1).....	25 Lbs

**TIRES/CASTORS**

Front Castor Assy.....	7"
Mid Wheel Tire Flat Free.....	10"
Rear Flat Free.....	7"

**BATTERIES**

BAT-U1 Module Capacity.....	12v35Ah x2
-----------------------------	------------

**CHARGER**

Charger Off Board.....	4 amp
------------------------	-------

**MOTOR POWER OUTPUT**

Horsepower.....	0.8 x2**
-----------------	----------

**DRIVING SPECS**

Turning Radius.....	15.75"
Range.....	up to 18.8 miles
Maximum Speed.....	5 mph
Maximum Climbing Angle.....	10% / 6°

All specifications are subject to change without notice. Shopprider Mobility Products, Inc® reserves the rights of any changes to the unit:  
 1) Include the anti-tip wheel or rear castor.  
 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.  
 3) The actual driving range varies with the factors shown below:  
 a) The weight of user.  
 b) Ground surface and terrain.  
 c) Battery capacity and conditions.  
 d) Type of charger.  
 e) Ambient Temperature.  
 f) The way and/or method of driving.  
 g) If the battery and mechanical moving parts are fully broken in.

**TOTAL \$**

\*\*STANDARD WITH UNIT