



## GK73 DASHER-9

Travel Mid-size  
Scooter  
with Splash n' Dash

### PRICE LIST & ORDER FORM Effective November 2006

For easy ordering and information contact Customer Service Toll Free at:

800-743-0772 Fax: 877-797-7081

Name _____	Phone No _____
Acct No _____ PO # _____	Date _____
Address _____	City/State _____ Zip _____
Ship To _____	Phone No _____
Address _____	City/State _____ Zip _____

**Stock** For your convenience, order by stock number (e.g. BAT-1/2)

#### BASE SCOOTER

- |                          |      |          |       |            |
|--------------------------|------|----------|-------|------------|
| <input type="checkbox"/> | GK73 | Dasher™9 | K0806 | \$2,195.00 |
|--------------------------|------|----------|-------|------------|
- \*Red, Blue, and Champagne Color Quick Change Panels Included*

#### BATTERY

- |                          |             |                                            |  |              |
|--------------------------|-------------|--------------------------------------------|--|--------------|
| <input type="checkbox"/> | Please Call | Battery Pack w/ 2 x 12v-21Ah SLA Batteries |  | TBA**        |
| <input type="checkbox"/> | BAT-1/2     | 12V - 21Ah SLA battery x2                  |  | \$90.00 ea** |
- Note: Complimentary battery has limited warranty

#### SEATS (W"xD"xH")

- |                          |              |                                          |  |            |
|--------------------------|--------------|------------------------------------------|--|------------|
| <input type="checkbox"/> | SF181614G    | 18"Wx16"D Stadium Style Fish-On Assy     |  | \$240.00** |
| <input type="checkbox"/> | 104215-77100 | Seat/Back Cushion Set-Vinyl Pads         |  | \$160.00** |
| <input type="checkbox"/> | SP161618RG   | Upgrade 16"w Seat Pan/Capt Back Recline* |  | \$150.00   |
| <input type="checkbox"/> | SP181821RG   | Upgrade 18"w Seat Pan/Capt Back Recline* |  | \$150.00   |

#### ARM REST ASSY

- |                          |                 |                                      |  |           |
|--------------------------|-----------------|--------------------------------------|--|-----------|
| <input type="checkbox"/> | 104300-69000-A3 | (L) Flip-up Arm Rest Assy w/ Arm Pad |  | \$56.00** |
| <input type="checkbox"/> | 104400-69000-A3 | (R) Flip-up Arm Rest Assy w/ Arm Pad |  | \$56.00** |

#### TIRES

- |                          |             |                                    |  |       |
|--------------------------|-------------|------------------------------------|--|-------|
| <input type="checkbox"/> | Please Call | Flat Free Tire Rear Wheel Assy x2  |  | TBA** |
| <input type="checkbox"/> | Please Call | Flat Free Tire Front Wheel Assy x1 |  | TBA** |

#### ANTI-TIP WHEELS

- |                          |              |                     |    |              |
|--------------------------|--------------|---------------------|----|--------------|
| <input type="checkbox"/> | 209928-89300 | Anti-Tip Wheel Assy | x2 | \$20.00 ea** |
|--------------------------|--------------|---------------------|----|--------------|

#### CHARGER

- |                          |              |                         |  |            |
|--------------------------|--------------|-------------------------|--|------------|
| <input type="checkbox"/> | 109201-88306 | 3 amp Charger Off-Board |  | \$180.00** |
|--------------------------|--------------|-------------------------|--|------------|

#### ACCESSORIES

- |                          |                 |                           |  |            |
|--------------------------|-----------------|---------------------------|--|------------|
| <input type="checkbox"/> | 209916-88100    | Wire Mesh Front Basket    |  | \$50.00**  |
| <input type="checkbox"/> | OXY-03          | Mesh Oxygen Tank Holder   |  | \$60.00    |
| <input type="checkbox"/> | CH-02-GK(K)     | Metal Oxygen Holder Kit   |  | \$150.00   |
| <input type="checkbox"/> | CCH-003-GK(K)   | Cane Holder Kit           |  | \$150.00   |
| <input type="checkbox"/> | RB-0125-GK(K)   | Rear basket Kit           |  | \$150.00   |
| <input type="checkbox"/> | WH-0185-GK(K)   | Walker Holder Kit         |  | \$150.00   |
| <input type="checkbox"/> | WIG-0168        | Wig-wag extensions (pair) |  | \$20.00    |
| <input type="checkbox"/> | 300909-03       | Safety Belt               |  | \$34.00    |
| <input type="checkbox"/> | 800904-69400-65 | GK73 Shroud in Champagne  |  | \$180.00** |
| <input type="checkbox"/> | 800904-69400-22 | GK73 Shroud in Red        |  | \$180.00** |
| <input type="checkbox"/> | 800904-69400-4A | GK73 Shroud in Blue       |  | \$180.00** |

#### DASHER 9 SPECIFICATIONS

##### OVERALL DIMENSIONS

Base Length.....	38.8"
Base Width.....	22"
Height.....	36"
(Ground to Seat Top)	
Ground Clearance.....	1.5"

##### WEIGHTS

User Weight Capacity.....	300 Lbs
Total Unit Weight.....	110 Lbs
Weight of Heaviest Part.....	35 Lbs
Battery Module Weight (Pack).....	35 Lbs

##### TIRES

Number, Size of Front Tire.....	1, 9"
Number, Size of Rear Tire.....	2, 9"

##### BATTERIES

BAT-12 Module Capacity.....	12v21Ah x2
-----------------------------	------------

##### CHARGER

Charger-Off Board.....	3 amp
------------------------	-------

##### MOTOR POWER OUTPUT

Horsepower.....	.75hp
-----------------	-------

##### DRIVING SPECS

Turning Radius.....	31.5"
Range*.....	7.5 miles
Maximum Speed.....	3.8 mph
Maximum Climbing Angle.....	10% (6°)

#### \*\*STANDARD WITH UNIT

#### TOTAL \$

All specifications are subject to change without notice.

Shoprider Mobility Products, Inc® reserves the rights of any changes to the unit:

- 1) Include the anti-tip wheel or rear castor.
- 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
- 3) \*The actual driving range varies with the factors shown below:
  - a) The weight of user.
  - b) Ground surface and terrain.
  - c) Battery capacity and conditions.
  - d) Type of charger.
  - e) Ambient Temperature.
  - f) The way and/or method of driving.
  - g) If the battery and mechanical moving parts are fully broken in.