

# PRICE LIST & ORDER FORM

# SHOPRIDER®

For easy ordering and information, please contact Customer Service **Toll Free at 1.800.743.0772 FAX: 1.877.797.7081** Effective July 2009



**GK73 DASHER 9**

NAME: \_\_\_\_\_  
 PHONE NO: \_\_\_\_\_  
 ACCT NO: \_\_\_\_\_  
 PO #: \_\_\_\_\_ DATE: \_\_\_\_\_  
 ADDRESS: \_\_\_\_\_  
 CITY: \_\_\_\_\_ STATE: \_\_\_\_\_ ZIP: \_\_\_\_\_  
 SHIP TO: \_\_\_\_\_  
 PHONE NO: \_\_\_\_\_  
 ADDRESS: \_\_\_\_\_  
 CITY: \_\_\_\_\_ STATE: \_\_\_\_\_ ZIP: \_\_\_\_\_

## SPECIFICATIONS

### OVERALL DIMENSIONS

Base Length.....38.5"  
 Base Width.....22.5"  
 Height.....36"  
 Ground to Seat Top  
 Ground Clearance.....1.5"

### WEIGHTS

User Weight Capacity.....300 Lbs  
 Total Unit Weight.....110 Lbs  
 Weight of Heaviest Part.....35 Lbs  
 Battery Weight Module (pack) .....35 Lbs

### TIRES

Number, Size of Front Tire.....1, 9"  
 Number, Size of Rear Tire.....2, 9"

### BATTERIES

BAT-21 Module Capacity.....12v21Ah x2

### CHARGER

Charger-Off Board.....2 amp

### MOTOR POWER OUTPUT

Horsepower.....0.75 hp

### DRIVING SPECS

Turning Radius .....31.5"  
 Range .....7.5 mph  
 Maximum Speed .....3.8 mph  
 Maximum Climbing Angle .....10% / 6°

### BASE SCOOTER

GK73

Dasher 9

K0806

### MSRP

\$2,195.00

\*Red, Blue, and Champagne Color Quick Change Panels Included

### BATTERY

BAT-21

12V - 21Ah SLA battery x2

\$110.00ea\*\*

Note: Complimentary battery has limited warranty

### SEATS (W"xD"xH")

SF181614G

18"Wx16"D Stadium Style Fish-On Assy

\$240.00\*\*

104215-77100

Seat/Back Cushion Set-Vinyl Pads

\$160.00\*\*

SP161618RG

Upgrade 16"w Seat Pan/Capt Back Recline\*

\$400.00

SP181821RG

Upgrade 18"w Seat Pan/Capt Back Recline\*

\$400.00

### ARM REST ASSY

104300-00106-A3

(L) Flipback High Adjustable Armrest

\$140.00\*\*

104400-00106-A3

(R) Flipback High Adjustable Armrest

\$140.00\*\*

### ANTI-TIP WHEELS

209928-89300

Anti-Tip Wheel Assy x2

\$20.00 ea\*\*

### CHARGER

109201-88803-A2

2 amp Charger Off-Board

\$100.00\*\*

### ACCESSORIES

209916-88100

Wire Mesh Front Basket

\$50.00\*\*

OXY-03

Mesh Oxygen Tank Holder

\$60.00

CH-02-GK(K)

Metal Oxygen Holder Kit

\$150.00

RB-0125-GK(K)

Rear basket Kit

\$150.00

WH-0185-GK(K)

Walker Holder Kit

\$150.00

WIG-0168

Wig-wag extensions (pair)

\$20.00

300909-03

Safety Belt

\$34.00

800904-69400-65

GK73 Shroud in Champagne

\$180.00\*\*

800904-69400-22

GK73 Shroud in Red

\$180.00\*\*

800904-69400-4A

GK73 Shroud in Blue

\$180.00\*\*

All specifications are subject to change without notice. Shoprider Mobility Products, Inc® reserves the rights of any changes to the unit:

- 1) Include the anti-tip wheel or rear castor.
- 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown below:
  - a) The weight of user.
  - b) Ground surface and terrain.
  - c) Battery capacity and conditions.
  - d) Type of charger.
  - e) Ambient Temperature.
  - f) The way and/or method of driving.
  - g) If the battery and mechanical moving parts are fully broken in.

\*\*STANDARD WITH UNIT

TOTAL \$ \_\_\_\_\_