

# PRICE LIST & ORDER FORM



For easy ordering and information, please contact Customer Service **Toll Free at 1.800.743.0772 FAX: 1.877.797.7081** Effective Jan 2011



**GK83-RED COOPER**

NAME: \_\_\_\_\_  
 PHONE NO: \_\_\_\_\_  
 ACCT NO: \_\_\_\_\_  
 PO #: \_\_\_\_\_ DATE: \_\_\_\_\_  
 ADDRESS: \_\_\_\_\_  
 CITY: \_\_\_\_\_ STATE: \_\_\_\_\_ ZIP: \_\_\_\_\_  
 SHIP TO: \_\_\_\_\_  
 PHONE NO: \_\_\_\_\_  
 ADDRESS: \_\_\_\_\_  
 CITY: \_\_\_\_\_ STATE: \_\_\_\_\_ ZIP: \_\_\_\_\_

## SPECIFICATIONS

### OVERALL DIMENSIONS

Base Length with Basket.....38"  
 Base Width.....22"  
 Height.....31-35"  
 Ground to Seat Top  
 Ground Clearance.....1.5"

### WEIGHTS

User Weight Capacity.....250 Lbs  
 Total Unit Weight.....87 Lbs  
 Weight of Heaviest Part.....25 Lbs  
 Battery Weight Module (pack) .....22 Lbs

### TIRES

Number, Size of Front Tire.....1, 8"  
 Number, Size of Rear Tire.....2, 8"

### BATTERIES

BAT-12 Module Capacity.....12v12Ah x2

### CHARGER

Charger-Off Board.....2 amp

### MOTOR POWER OUTPUT

Horsepower.....0.5 hp

### DRIVING SPECS

Turning Radius .....33"  
 Range .....7.5 mph  
 Maximum Speed .....3.8 mph  
 Maximum Climbing Angle .....10% / 6°

All specifications are subject to change without notice. Shoprider Mobility Products, Inc.® reserves the rights of any changes to the unit:

- 1) Include the anti-tip wheel or rear castor.
- 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown below:
  - a) The weight of user.
  - b) Ground surface and terrain.
  - c) Battery capacity and conditions.
  - d) Type of charger.
  - e) Ambient Temperature.
  - f) The way and/or method of driving.
  - g) If the battery and mechanical moving parts are fully broken in.

### BASE SCOOTER

GK83-RED Cooper MSRP \$1,795.00

### BATTERY

800501-695-01 Battery Pack w/ 2 x 12v-12amp SLA Batteries \$300.00\*\*  
 BAT-12 12V - 12amp SLA battery x2 \$75.00 ea\*\*

### SEATS (W"xD"xH")

SF181614G(K) 18"Wx16"D Stadium Style Fish-On Assy \$160.00\*\*  
 104213-66700-A2 Seat/Back Cushion Set-Vinyl Pads \$80.00\*\*  
 SP161618RG Upgrade 16"w Seat Pan/Capt Back Recline\* \$400.00  
 \*Also avail in 18"w. Seat cushion not sold.

### ARM REST ASSY

104300-69000-A3 (L) Flip-up Arm Rest Assy w/ Arm Pad \$56.00\*\*  
 104400-69000-A3 (R) Flip-up Arm Rest Assy w/ Arm Pad \$56.00\*\*

### TIRES

105201-69101 Solid Tire Rear Wheel Assy x2 \$90.00 ea\*\*  
 105101-69102 Solid Tire Front Wheel Assy x1 \$100.00ea\*\*

### ANTI-TIP WHEELS

209928-89300 Anti-Tip Wheel Assy x2 \$20.00 ea\*\*

### CHARGER

109201-88800-A3 2 amp Charger Off-Board \$100.00\*\*

### ACCESSORIES

209916-78700 Wire Mesh Front Basket \$40.00\*\*  
 OXY-03 Mesh Oxygen Tank Holder \$60.00  
 0210(K) Metal Oxygen Holder Kit \$150.00  
 CCH-003-GK(K) Cane Holder Kit \$150.00  
 WH-0185-GK(K) Walker Holder Kit \$150.00  
 WIG-0168 Wig-wag extensions (pair) \$20.00  
 300909-03 Safety Belt \$34.00  
 800904-69501-4A Shroud covers in blue \$140.00  
 800904-69501-65 Shroud covers in champagne \$140.00  
 800904-69501-13 Shroud covers in dark grey \$140.00

\*\*STANDARD WITH UNIT

TOTAL \$ \_\_\_\_\_