



GK83
DASHER-3
 Travel Scooter
 with Splash n' Dash

PRICE LIST & ORDER FORM Effective November 2006

For easy ordering and information contact Customer Service Toll Free at:

800-743-0772 Fax: 877-797-7081

| | |
|---------------|----------------------------|
| Name _____ | Phone No _____ |
| Acct No _____ | PO # _____ |
| Date _____ | |
| Address _____ | City/State _____ Zip _____ |
| Ship To _____ | Phone No _____ |
| Address _____ | City/State _____ Zip _____ |

Stock For your convenience, order by stock number (e.g. BAT-12)

BASE SCOOTER

GK83 Dasher™3 K0800 \$1,795.00

*Red, Blue, and Champagne Color Quick Change Panels Included

BATTERY

800501-695-01 Battery Pack w/ 2 x 12v-12amp SLA Batteries \$300.00**
 BAT-12 12V - 12amp SLA battery x2 \$75.00 ea**

*Complimentary Battery not warranted
 18"Wx16"D Stadium Style Fish-On Assy \$160.00**
 Seat/Back Cushion Set-Vinyl Pads \$80.00**
 Upgrade 16"w Seat Pan/Capt Back Recline* \$400.00

SEATS (W"xD"xH")

SF181614G(K) 18"Wx16"D Stadium Style Fish-On Assy \$160.00**
 104213-66700-A2 Seat/Back Cushion Set-Vinyl Pads \$80.00**
 SP161618RG Upgrade 16"w Seat Pan/Capt Back Recline* \$400.00

ARM REST ASSY

104300-69000-A3 (L) Flip-up Arm Rest Assy w/ Arm Pad \$56.00**
 104400-69000-A3 (R) Flip-up Arm Rest Assy w/ Arm Pad \$56.00**

*Also avail in 18"w. Seat cushion not sold.

TIRES

105201-69101 Solid Tire Rear Wheel Assy x2 \$90.00 ea**
 105101-69102 Solid Tire Front Wheel Assy x1 \$100.00 ea**

ANTI-TIP WHEELS

209928-89300 Anti-Tip Wheel Assy x2 \$20.00 ea**

CHARGER

109201-88800-A3 2 amp Charger Off-Board \$100.00**

ACCESSORIES

209916-78700 Wire Mesh Front Basket \$40.00**
 OXY-03 Mesh Oxygen Tank Holder \$60.00
 CH-02-GK(K) Metal Oxygen Holder Kit \$150.00
 CCH-003-GK(K) Cane Holder Kit \$150.00
 WH-0185-GK(K) Walker Holder Kit \$150.00
 WIG-0168 Wig-wag extensions (pair) \$20.00
 300909-03 Safety Belt \$34.00
 CVR-GK83-CAM Camouflage Color Changing Panel \$300.00
 CVR-GK83-FLAG Stars & Strips Color Changing Panel \$300.00
 CVR-GK83-FLAME Flames Color Changing Panel \$300.00

DASHER 3 SPECIFICATIONS

OVERALL DIMENSIONS

Base Length..... 42.5"
 Base Width..... 22"
 Height..... 35"
 (Ground to Seat Top)
 Ground Clearance..... 1.5"

WEIGHTS

User Weight Capacity..... 250 Lbs
 Total Unit Weight..... 89 Lbs
 Weight of Heaviest Part..... 26 Lbs
 Battery Module Weight (Pack)..... 22 Lbs

TIRES

Number, Size of Front Tire..... 1, 8"
 Number, Size of Rear Tire..... 2, 8"

BATTERIES

BAT-12 Module Capacity..... 12v12Ah x2

CHARGER

Charger-Off Board..... 2 amp

MOTOR POWER OUTPUT

Horsepower..... .5hp

DRIVING SPECS

Turning Radius..... 33"
 Range..... 7.5 miles
 Maximum Speed..... 3.8 mph
 Maximum Climbing Angle..... 10% (6°)

****STANDARD WITH UNIT**

TOTAL \$

All specifications are subject to change without notice.

Shopperider Mobility Products, Inc® reserves the rights of any changes to the unit:

- 1) Include the anti-tip wheel or rear castor.
- 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown below:
 - a) The weight of user.
 - b) Ground surface and terrain.
 - c) Battery capacity and conditions.
 - d) Type of charger.
 - e) Ambient Temperature.
 - f) The way and/or method of driving.
 - g) If the battery and mechanical moving parts are fully broken in.