



**GK8
DASHER-4**

**Travel Scooter
with Splash n' Dash**

PRICE LIST & ORDER FORM Effective November 2006

For easy ordering and information contact Customer Service Toll Free at:
800-743-0772 Fax: 877-797-7081

Name _____	Phone No _____
Acct No _____	PO # _____
Address _____	Date _____
Ship To _____	City/State _____ Zip _____
Address _____	Phone No _____
City/State _____	Zip _____

Stock For your convenience, order by stock number (e.g. BAT-12)

BASE SCOOTER

<input type="checkbox"/>	GK8	Dasher™4	K0800	\$1,895.00
--------------------------	-----	----------	-------	------------

*Red, Blue, and Champagne Color Quick Change Panels Included

BATTERY

<input type="checkbox"/>	800501-695-01	Battery Pack w/ 2 x 12v-12amp SLA Batteries		\$300.00**
<input type="checkbox"/>	BAT-12	12V - 12amp SLA battery	x2	\$70.00 ea**

SEATS (W"xD"xH")

<input type="checkbox"/>	SF181614G(K)	18"Wx16"D Stadium Style Fish-On Assy		\$160.00**
<input type="checkbox"/>	104213-66700-A2	Seat/Back Cushion Set-Vinyl Pads		\$80.00**
<input type="checkbox"/>	SP161618RG	Upgrade 16"w Seat Pan/Capt Back Recline*		\$150.00

*Complimentary Battery not warranted
*Also avail in 18"w. Seat cushion not sold.

ARM REST ASSY

<input type="checkbox"/>	104300-69100-A3	(L) Flip-up Arm Rest Assy w/ Arm Pad		\$56.00**
<input type="checkbox"/>	104300-69100-A4	(R) Flip-up Arm Rest Assy w/ Arm Pad		\$56.00**

TIRES

<input type="checkbox"/>	105201-69101	Solid Tire Rear Wheel Assy	x2	\$90.00 ea**
<input type="checkbox"/>	105101-69501	Solid Tire Front Wheel Assy	x2	\$90.00 ea**

ANTI-TIP WHEELS

<input type="checkbox"/>	209928-89300	Anti-Tip Wheel Assy	x2	\$20.00 ea**
--------------------------	--------------	---------------------	----	--------------

CHARGER

<input type="checkbox"/>	109201-88805	2 amp Charger Off-Board		\$100.00**
--------------------------	--------------	-------------------------	--	------------

ACCESSORIES

<input type="checkbox"/>	209916-78700	Wire Mesh Front Basket		\$40.00**
<input type="checkbox"/>	OXY-03	Mesh Oxygen Tank Holder		\$60.00
<input type="checkbox"/>	CH-02-GK(K)	Metal Oxygen Holder Kit		\$150.00
<input type="checkbox"/>	CCH-003-GK(K)	Cane Holder Kit		\$150.00
<input type="checkbox"/>	RB-0125-GK(K)	Rear basket Kit		\$150.00
<input type="checkbox"/>	WH-0185-GK(K)	Walker Holder Kit		\$150.00
<input type="checkbox"/>	WIG-0168	Wig-wag extensions (pair)		\$20.00
<input type="checkbox"/>	300909-03	Safety Belt		\$34.00
<input type="checkbox"/>	CVR-GK8-CAM	Camouflage Color Changing Panel		\$300.00
<input type="checkbox"/>	CVR-GK8-FLAG	Stars & Strips Color Changing Panel		\$300.00
<input type="checkbox"/>	CVR-GK8-FLAME	Flames Color Changing Panel		\$300.00

DASHER 4 SPECIFICATIONS

OVERALL DIMENSIONS

Base Length.....	44"
Base Width.....	22"
Height.....	33"
(Ground to Seat Top)	
Ground Clearance.....	1.5"

WEIGHTS

User Weight Capacity.....	250 Lbs
Total Unit Weight.....	94 Lbs
Weight of Heaviest Part.....	33 Lbs
Battery Module Weight (Pack).....	22 Lbs

TIRES

Number, Size of Front Tire.....	2, 7"
Number, Size of Rear Tire.....	2, 8"

BATTERIES

BAT-12 Module Capacity.....	12v12Ah x2
-----------------------------	------------

CHARGER

Charger-Off Board.....	2 amp
------------------------	-------

MOTOR POWER OUTPUT

Horsepower.....	.5
-----------------	----

DRIVING SPECS

Turning Radius.....	37"
Range.....	7.5 miles
Maximum Speed.....	3.8 mph
Maximum Climbing Angle.....	10% (6°)

TOTAL \$

All specifications are subject to change without notice.

Shoprider Mobility Products, Inc® reserves the rights of any changes to the unit:

- 1) Include the anti-tip wheel or rear castor.
- 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown below:
 - a) The weight of user.
 - b) Ground surface and terrain.
 - c) Battery capacity and conditions.
 - d) Type of charger.
 - e) Ambient Temperature.
 - f) The way and/or method of driving.
 - g) If the battery and mechanical moving parts are fully broken in.