



2/13/09

PRICE LIST & ORDER FORM Effective November 2006

For easy ordering and information contact Customer Service Toll Free at:

800-743-0772 Fax: 877-797-7081

Name _____	Phone No _____
Acct No _____	PO # _____
Address _____	Date _____
City/State _____	Zip _____
Ship To _____	Phone No _____
Address _____	City/State _____
	Zip _____

Stock For your convenience, order by stock number (e.g. BAT-12)

BASE SCOOTER

TE-SL7-3 ECHO 3® RED/BLUE \$1,495.00

BATTERY

800501-699-00 Battery Pack w/ 2 x 12v-10amp SLA Batteries \$280.00**
 BAT-10 12V - 10amp SLA battery x2 \$70.00 ea**
 *Complimentary Battery not warranted

SEATS (W"xD"xH")

SF171413B 17"Wx14"D Stadium Style Fish-On w/PAD PLEASE CALL**
 ECHO-UP(K) Upgrade KIT w/armrest & back pad pouch \$80.00
 SF181614G(K) Upgrade 18"Wx16"D Fish on w/Renoveable Pads \$160.00
 SP161618RG Upgrade 16"Wx16"D Pan Seat 2-tone grey recline \$400.00

ARM REST ASSY

ECHO-UP(K) Upgrade KIT w/armrest & back pad pouch \$80.00

TIRES

105201-69101 Solid Tire Rear Wheel Assy x2 \$90.00 ea**
 105101-69501 Solid Tire Front Wheel Assy x1 \$90.00 ea**

ANTI-TIP WHEELS

209928-89300 Anti-Tip Wheel Assy x2 \$20.00 ea**

CHARGER

109201-88800-A3 2 amp Charger Off-Board \$100.00**

ACCESSORIES

209916-78700 Wire Mesh Front Basket \$40.00**
 OXY-03 Mesh Oxygen Tank Holder \$60.00
 WIG-0168 Wig-wag extensions (pair) \$20.00
 300909-03 Safety Belt \$34.00

TOTAL \$

**STANDARD WITH UNIT

All specifications are subject to change without notice.

Shoprider Mobility Products, Inc® reserves the rights of any changes to the unit:

- 1) Include the anti-tip wheel or rear castor.
- 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown below:
 - a) The weight of user.
 - b) Ground surface and terrain.
 - c) Battery capacity and conditions.
 - d) Type of charger.
 - e) Ambient Temperature.
 - f) The way and/or method of driving.
 - g) If the battery and mechanical moving parts are fully broken in.

ECHO 3 SPECIFICATIONS

OVERALL DIMENSIONS

Base Length..... 37"
 Base Width..... 20"
 Height..... 36"
 (Ground to Seat Top)
 Ground Clearance..... 1.5"

WEIGHTS

User Weight Capacity..... 250 Lbs
 Total Unit Weight..... 70 Lbs
 Weight of Heaviest Part..... 35 Lbs
 Battery Module Weight (Pack)... 15 Lbs

TIRES

Number, Size of Front Tire..... 1, 8"
 Number, Size of Rear Tire..... 2, 8"

BATTERIES

BAT-10 Module Capacity..... 12v10Ah x2

CHARGER

Charger-Off Board..... 2 amp

MOTOR POWER OUTPUT

Horsepower..... .5

DRIVING SPECS

Turning Radius..... 32"
 Range..... 8 miles
 Maximum Speed..... 3.8 mph
 Maximum Climbing Angle..... 10% (6°)