

PRICE LIST & ORDER FORM



For easy ordering and information, please contact Customer Service **Toll Free at 1.800.743.0772 FAX: 1.877.797.7081** Effective Oct,2007



889XLSN FLAGSHIP

NAME: _____
 PHONE NO: _____
 ACCT NO: _____
 PO #: _____ DATE: _____
 ADDRESS: _____
 CITY: _____ STATE: _____ ZIP: _____
 SHIP TO: _____
 PHONE NO: _____
 ADDRESS: _____
 CITY: _____ STATE: _____ ZIP: _____

SPECIFICATIONS

OVERALL DIMENSIONS

Base Length.....63"
 Base Width.....29"
 Height.....64"
 Ground to Seat Top
 Ground Clearance.....5"

STANDARD SEAT

Black vinyl, semi-recline, high back, adjustable height head rest
 Seat W x D x H.....20"x18"x20"
 Seat Pad Cushion Height.....4"

WEIGHTS

User Weight Capacity.....350 Lbs
 Total Unit Weight.....328 Lbs
 Battery Module Weight.....114 Lbs

TIRES

Front Wheel.....12"
 Rear Flat Free.....12"

BATTERIES

BAT-75 Module Capacity.....12v75Ah x2

CHARGER

Charger-Off Board.....8 amp

MOTOR POWER OUTPUT

Horsepower.....2.0 hp

DRIVING SPECS

Turning Radius70"
 Rangeup to 25 miles
 Maximum Speed7 mph
 Maximum Climbing Angle15% / 9°

BASE SCOOTER

889XLSN Flagship Scooter
 RED/BLUE

MSRP

\$8,995.00

BATTERY

BAT-75F 12V 75AMP x2 \$600.00

SEATS (W"xD"xH")

SV201820RB 20"x18" Pillow-Top Recline Seat-Vinyl \$370.00**
 No upgrades seat available for this model

TIRES

TAR330100 Rear 12" Pnuematic Tires-889B-4,778XLSBN,889XLSBN \$50.00**
 TAF330100 Front 12" Pnuematic Tires-889B-4,778XLSBN,889XLSBN \$50.00**
 205102-89100 Front & Rear 12" Tire Tubes-889B-4,778XLSBN,889XLSBN \$20.00**
 SF-330100-SC(K) Front Wheel assy kit; 889B-4 solid front \$248.00
 SR-330100XL4(K) Rear Wheel assy kit; 889B-4 solid rear \$248.00

CHARGER

109201-88500 8 amp Charger Off-Board (889X,778XLSBN,889XLSBN) \$300.00**

ACCESSORIES

300909-03 Safety Belt \$34.00
 5008-10020-212 Tightening knob for accessory X2 \$20.00
 202908-88105 Flagship side mirror \$18.00
 800920-67200-10 Rear basket assembly w/o cover \$250.00
 800920-6720-10 Lockable box assembly w/cover \$300.00

**STANDARD WITH UNIT

TOTAL \$ _____

All specifications are subject to change without notice.
 Shoprider Mobility Products, Inc® reserves the rights of any changes to the unit:

- 1) Include the anti-tip wheel or rear castor.
- 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown below:
 - a) The weight of user.
 - b) Ground surface and terrain.
 - c) Battery capacity and conditions.
 - d) Type of charger.
 - e) Ambient Temperature.
 - f) The way and/or method of driving.
 - g) If the battery and mechanical moving parts are fully broken in.