

PRICE LIST & ORDER FORM



For easy ordering and information, please contact Customer Service **Toll Free at 1.800.743.0772 FAX: 1.877.797.7081** Effective July 2009



FPC11 FOLDING CHAIR

NAME: _____
 PHONE NO: _____
 ACCT NO: _____
 PO #: _____ DATE: _____
 ADDRESS: _____
 CITY: _____ STATE: _____ ZIP: _____
 SHIP TO: _____
 PHONE NO: _____
 ADDRESS: _____
 CITY: _____ STATE: _____ ZIP: _____

SPECIFICATIONS

OVERALL DIMENSIONS [18" , 20"]

Base Length.....46"
 Base Width.....25/27"
 Height.....37"
 (Ground to Seat Top)
 Ground Clearance.....3"

WEIGHTS

User Weight Capacity.....300 Lbs
 Total Unit Weight w/Batteries.....130 Lbs
 Weight of the Heaviest Part.....80 Lbs
 Battery Module Weight (U1).....25 Lbs x2

TIRES/CASTORS

Front Castor Assy.....8"
 Rear Tire-Flat Free.....12"

BATTERIES

BAT-U1 Module Capacity.....12V35Ahx2

CHARGER

Charger-Off Board.....5 amp

MOTOR POWER OUTPUT

Horsepower.....0.60 x2

DRIVING SPECS

Turning Radius24"
 Range10-12.5 miles
 Maximum Speed5 mph
 Maximum Climbing Angle14% / 8°

All specifications are subject to change without notice.
 Shoprider Mobility Products, Inc.® reserves the rights of any changes to the unit:

- 1) Include the anti-tip wheel or rear castor.
- 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown below:
 - a) The weight of user.
 - b) Ground surface and terrain.
 - c) Battery capacity and conditions.
 - d) Type of charger.
 - e) Ambient Temperature.
 - f) The way and/or method of driving.
 - g) If the battery and mechanical moving parts are fully broken in.

BASE POWERCHAIR

PHFW-1118 FPC-1118 BGRD 18" wide K0822 \$4,895.00
 PHFW-1120 FPC-1120 BGRD 20" wide K0822 \$4,895.00

BATTERY

BAT-U1 12V 35AMP x2 \$140.00**ea

BATTERY ACCESSORIES

800501-007 Complete Battery Kit x1 \$425.00**

SEAT ACCESSORIES

204201-00700 18" Seat Cushion \$44.00
 204201-00701 20" Seat Cushion \$44.00
 SR263861 18 x 16 General Seat Cushion E2601 \$79.00
 SR263081 20 x 18 General Seat Cushion E2601 \$79.00

ARM REST ASSY

104300-00201 (L) Upper Arm Rest Assy \$32.00**
 1319-00100-00 (L) Lower Arm Rest Assy \$32.00**
 104400-00201 (R) Upper Arm Rest Assy \$32.00**
 1319-00100-01 (R) Lower Arm Rest Assy \$32.00**

DRIVE TIRES (Rear Wheel)

105201-00700-50 12" Flat-Free Wheel Assy x2 \$200.00**ea

CASTORS (Front Wheel)

1315-00202-00 Front Fork Assy x2 \$120.00**ea
 105101-00701 8" Flat-Free Front Wheel Assy x2 \$140.00**ea

ANTI-TIP WHEELS

209928-00700 Anti-Tip Wheel Assy x2 \$20.00**

FOOT REST

ELRFPC-R(K) Elevating Leg Rest Kit (L&R) \$160.00
 ELRFPC-B(K) Elevating Leg Rest Kit Black (L&R) \$160.00
 1329-00404-00 Swing Away Foot Rest (R) \$80.00**
 1329-00403-00 Swing Away Foot Rest (L) \$80.00**

CONTROLLER

1361-00206-00 VSI Joystick Controller \$998.00**
 1343-PHFW-SQ VSI Controller Bracket \$36.00**

CHARGER

109201-88124 5 amp Charger Off-Board \$230.00**

ACCESSORIES

OXY-03 Mesh Oxygen Tank Holder \$60.00
 300909-03 Safety Belt \$34.00**

**STANDARD WITH UNIT

TOTAL \$ _____