



**889B-4**  
Sprinter XL4

**889-3XLS**  
Sprinter Jumbo 3

**889-4XLS**  
Sprinter XL4 Deluxe

**889X**  
Sprinter Jumbo XL

**PRICE LIST & ORDER FORM** Effective November 2006

For easy ordering and information  
contact Customer Service Toll Free at:

**800-743-0772 Fax: 877-797-7081**

Name _____	Phone No _____
Acct No _____	PO # _____
Date _____	
Address _____	City/State _____
Zip _____	
Ship To _____	Phone No _____
Address _____	City/State _____
Zip _____	

**BASE SCOOTER**

<input type="checkbox"/>	889B-4	SPRINTER XL4	Burgundy	K0899	\$2,595.00
<input type="checkbox"/>	778XLSBN	SPRINTER XL3 DELUXE	Burgundy	K0899	\$5,295.00
<input type="checkbox"/>	889XLSBN	SPRINTER XL4 DELUXE	Burgundy	K0899	\$5,995.00
<input type="checkbox"/>	889X	SPRINTER JUMBO XL	Burgundy	K0899	\$5,995.00

**BATTERY**

<input type="checkbox"/>	BAT-U1	x2	12V 35AMP-889B-4	\$220.00**
<input type="checkbox"/>	BAT-55	X2	12V 55AMP-778XLSBN,-889XLSBN,889X	\$380.00
<input type="checkbox"/>	BAT-75	x2	12V 75AMP-778XLSBN,-889XLSBN,889X	\$500.00

COMPLIMENTARY BATTERIES ARE NOT WARRANTED

**SEATS (W"xD"xH")**

<input type="checkbox"/>	SV171/21815G	17½x18" Deluxe Captain Seat Vinyl-889B-4	\$408.00**
<input type="checkbox"/>	SC171/21814G	Upgrade 17½" Deluxe Capt Seat Grey Cloth	\$150.00
<input type="checkbox"/>	SV171/21714RB	Upgrade 17½" Pillow-Top Recline Seat Vinyl	\$150.00
<input type="checkbox"/>	SV201820RG/B	20"x18" Highback ReclineSeat-Vinyl-778XLSBN,889XLSBN	\$370.00**
<input type="checkbox"/>	SV201820RB	Upgrade 20"x18" Pillow-Top Recline Seat-Vinyl	\$150.00
<input type="checkbox"/>	SC201822RG	Upgrade 20"x18" Recline Highback Seat-Cloth	\$150.00
<input type="checkbox"/>	SC221922RG	Upgrade 22"x19" Highback Captain Seat-Cloth	\$150.00
<input type="checkbox"/>	SV232020B	23" Pillow-Top Black Seat-Vinyl-889X	\$450.00**

**ARM REST ASSY**

<input type="checkbox"/>	104300-88800	(L) Arm Rest Assy w/ Arm Pad	\$160.00**
<input type="checkbox"/>	104400-88800	(R) Arm Rest Assy w/ Arm Pad	\$160.00**

**TIRES**

<input type="checkbox"/>	TAR330100	Rear	12" Pnuematic Tires-889B-4,778XLSBN,889XLSBN	\$50.00**
<input type="checkbox"/>	TAF330100	Front	12" Pnuematic Tires-889B-4,778XLSBN,889XLSBN	\$50.00**
<input type="checkbox"/>	205102-89100	Front & Rear	12" Tire Tubes-889B-4,778XLSBN,889XLSBN	\$20.00**
<input type="checkbox"/>	STU9X (K)	Front & Rear	14" Foam Filled Tires GREY-889X	\$320.00
<input type="checkbox"/>	STU9DX(K)	Front & Rear	Solid Tire Upgrade-889B-4, 889XLSBN	\$320.00
<input type="checkbox"/>	STU-3XLS(K)	Front & Rear	Solid Tire Upgrade-778XLSBN	\$240.00

**CHARGER**

<input type="checkbox"/>	109201-88306	3 amp Charger Off-Board (889B-4)	\$180.00**
<input type="checkbox"/>	109201-88500	8 amp Charger Off-Board (889X,778XLSBN,889XLSBN)	\$300.00**

**ACCESSORIES**

<input type="checkbox"/>	209916-88100	Front Wire Basket (Large) -889X	\$50.00**
<input type="checkbox"/>	109901-88101	Front Wicker Basket-889B-4,778XLSBN,889XLSBN	\$70.00**
<input type="checkbox"/>	204506-67200	Accessory Adaptor ONLY REQUIRED-889X	\$140.00
<input type="checkbox"/>	CCH-003(K)	Cane Holder	\$150.00
<input type="checkbox"/>	CH-02(K)	Oxygen Tank Holder (Metal)	\$150.00
<input type="checkbox"/>	WH-0185(K)	Walker Holder	\$150.00
<input type="checkbox"/>	RB-0125(K)	Rear Wire Mesh Basket 21"Wx15"Lx9"D	\$170.00
<input type="checkbox"/>	300909-03	Safety Belt	\$34.00

**\*\*STANDARD WITH UNIT**

**TOTAL \$**

All specifications are subject to change without notice.  
Shoprider Mobility Products, Inc® reserves the rights of any changes to the unit:  
1) Include the anti-tip wheel or rear castor.  
2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.  
3) The actual driving range varies with the factors shown below:  
a) The weight of user.  
b) Ground surface and terrain.  
c) Battery capacity and conditions.  
d) Type of charger.  
e) Ambient Temperature.  
f) The way and/or method of driving.  
g) If the battery and mechanical moving parts are fully broken in.

**HEAVY DUTY SPECIFICATIONS**

**SPRINTER XL4**

Base LxWxH".....	52"x23"x30"
*Ground to Seat Top	
Ground Clearance.....	3"
User Weight Capacity.....	350 Lbs
Total Unit Weight.....	181 Lbs
Battery Module Weight (40Ah).....	32 Lbs x2
Number, Size Tire (front/rear).....	2,12"/2,12"
BAT-40 Module Capacity.....	12v40Ah x2
Charger-Off Board.....	3 amp
Turning Radius.....	54"
Range.....	up to 25 miles
Maximum Speed.....	5 mph
Maximum Climbing Angle.....	14% / 8°

**SPRINTER XL3 DELUXE**

Base LxWxH".....	51"x27"x35"
*Ground to Seat Top	
Ground Clearance.....	4"
User Weight Capacity.....	500 Lbs
Total Unit Weight.....	222 Lbs
Battery Module Weight (BAT-75).....	50 Lbs x2
Number, Size Tire (front/rear).....	1,12"/2,12"
BAT-75 Module Capacity.....	12v75Ah x2
Charger-Off Board.....	8 amp
Turning Radius.....	40"
Range.....	up to 25 miles
Maximum Speed.....	5 mph
Maximum Climbing Angle.....	14% / 8°

**SPRINTER XL4 DELUXE**

Base LxWxH".....	53"x27"x34"
*Ground to Seat Top	
Ground Clearance.....	4"
User Weight Capacity.....	500 Lbs
Total Unit Weight.....	240 Lbs
Battery Module Weight (BAT-75).....	50 Lbs x2
Number, Size Tire (front/rear).....	2,12"/2,12"
BAT-75 Module Capacity.....	12v75Ah x2
Charger-Off Board.....	8 amp
Turning Radius.....	72"
Range.....	up to 25 mile
Maximum Speed.....	5 mph
Maximum Climbing Angle.....	14% / 8°

**SPRINTER JUMBO XL**

Base LxWxH".....	56"x25"x37"
*Ground to Seat Top	
Ground Clearance.....	5"
User Weight Capacity.....	500 Lbs
Total Unit Weight.....	154 Lbs
Battery Module Weight (BAT-75).....	50 Lbs x2
Number, Size Tire (front/rear).....	2,14"/2,14"
BAT-75 Module Capacity.....	12v75Ah x
Charger-Off Board.....	8 amp
Turning Radius.....	72"
Range.....	up to 49 r
Maximum Speed.....	8 mph
Maximum Climbing Angle.....	14% / 8°