



FPC-11
Folding
Powerchair

PRICE LIST & ORDER FORM Effective November 2006
For easy ordering and information contact Customer Service Toll Free at:
800-743-0772 Fax: 877-797-7081

Name _____	Phone No _____
Acct No _____	PO # _____ Date _____
Address _____	City/State _____ Zip _____
Ship To _____	Phone No _____
Address _____	City/State _____ Zip _____

Stock For your convenience, order by stock number (e.g. BAT-U1)

BASE CHAIR				
<input type="checkbox"/>	PHFW-1118	FPC-1118 BGRD 18" wide	K0822	\$4,895.00
<input type="checkbox"/>	PHFW-1120	FPC-1120 BGRD 20" wide	K0822	\$4,895.00
BATTERY				
<input type="checkbox"/>	BAT-U1	12V 35AMP	x2	\$220.00
BATTERY ACCESSORIES				
<input type="checkbox"/>	209922-00700	Battery Lid w/ Cables	x2	\$72.00**
<input type="checkbox"/>	8306-00100-00	18" Battery Basket	x2	\$140.00**
<input type="checkbox"/>	8306-00100-01	20" Battery Basket	x2	\$140.00**
SEAT ACCESSORY				
<input type="checkbox"/>	204201-00700	18" Seat Cushion		\$44.00**
<input type="checkbox"/>	204201-00701	20" Seat Cushion		\$44.00**
ARM REST ASSY				
<input type="checkbox"/>	104300-00200	(L) Upper Arm Rest Assy		\$32.00**
<input type="checkbox"/>	1319-00100-00	(L) Lower Arm Rest Assy		\$32.00**
<input type="checkbox"/>	104400-00200	(R) Upper Arm Rest Assy		\$32.00**
<input type="checkbox"/>	1319-00100-01	(R) Lower Arm Rest Assy		\$32.00**
DRIVE TIRES (Rear Wheel)				
<input type="checkbox"/>	105201-00700-50	12" Flat-Free Wheel Assy	x2	\$400.00**
CASTORS (Front Wheel)				
<input type="checkbox"/>	232901-PHFW	Front Fork Assy	x2	\$88.00**
<input type="checkbox"/>	105101-00701	8" Flat-Free Front Wheel Ass x2		\$280.00**
ANTI-TIP WHEELS				
<input type="checkbox"/>	1329-00104-00	Anti-Tip Wheel Assy	x2	\$96.00**
FOOT REST				
<input type="checkbox"/>	ELRFPC-R(K)	Elevating Leg Rest Kit (L&R)		\$152.00
<input type="checkbox"/>	1329-00404-00	Swing Away Foot Rest (R)		\$80.00**
<input type="checkbox"/>	1329-00403-00	Swing Away Foot Rest (L)		\$80.00**
CONTROLLER				
<input type="checkbox"/>	1361-00206-00	VSI Joystick Controller		\$998.00**
<input type="checkbox"/>	1343-PHFW-SQ	VSI Controller Bracket		\$36.00**
CHARGER				
<input type="checkbox"/>	109201-88109-A2	4 amp Charger Off-Board		\$200.00**
ACCESSORIES				
<input type="checkbox"/>	OXY-03	Mesh Oxygen Tank Holder		\$60.00
<input type="checkbox"/>	ARM-BAG LEFT	Arm Pouch (left)		\$20.00
<input type="checkbox"/>	ARM-BAG RIGHT	Arm Pouch (right)		\$20.00
<input type="checkbox"/>	300909-03	Safety Belt		\$34.00**

FPC SPECIFICATIONS	
OVERALL DIMENSIONS	
Base Length.....	43"
Base Width.....	24"/26"
Height.....	37"
(Ground to Seat Top)	
Ground Clearance.....	3"
WEIGHTS	
User Weight Capacity.....	300 Lbs
Total Unit Weight w/ Batteries.....	130 Lbs
Weight of Heaviest Part.....	80 Lbs
Battery Module Weight (U1).....	25 Lbs x 2
TIRES/CASTORS	
Front Castor Assy.....	8"
Rear Tire-Flat Free.....	12"
BATTERIES	
BAT-U1 Module Capacity.....	12v35Ah x2
CHARGER	
Charger-Off Board.....	4 amp
MOTOR POWER OUTPUT	
Horsepower.....	0.60 x2
DRIVING SPECS	
Turning Radius.....	24"
Range.....	up to 20 miles
Maximum Speed.....	5 mph
Maximum Climbing Angle.....	14% / 8°

All specifications are subject to change without notice. Shopprider Mobility Products, Inc® reserves the rights of any changes to the unit:
 1) Include the anti-tip wheel or rear castor.
 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
 3) The actual driving range varies with the factors shown below:
 a) The weight of user.
 b) Ground surface and terrain.
 c) Battery capacity and conditions.
 d) Type of charger.
 e) Ambient Temperature.
 f) The way and/or method of driving.
 g) If the battery and mechanical moving parts are fully broken in.

TOTAL \$

**STANDARD WITH UNIT