



TE-787NA
SCOOTIE

PRICE LIST & ORDER FORM Effective November 2006

For easy ordering and information contact Customer Service Toll Free at:

800-743-0772 Fax: 877-797-7081

Name _____	Phone No _____
Acct No _____ PO # _____	Date _____
Address _____ City/State _____ Zip _____	
Ship To _____	Phone No _____
Address _____ City/State _____ Zip _____	

Stock For your convenience, order by stock number (e.g. BAT-12)

BASE SCOOTER

<input type="checkbox"/>	TE-787NA	SCOOTIE	RED	K0800	\$1,595.00
<input type="checkbox"/>	TE-787NA	SCOOTIE	BLUE	K0800	\$1,595.00

BATTERY

<input type="checkbox"/>	BAT-12	12V 12AMP	x2	\$140.00**
--------------------------	--------	-----------	----	------------

SEATS (W"xD"xH")

<input type="checkbox"/>	SF151413G	Grey Vinyl Fish-On Seat		\$100.00**
<input type="checkbox"/>	104101-78700	ABS Plastic Seat w/ Seat Pad		\$75.00

ARM REST ASSY

<input type="checkbox"/>	104300-78700	(L) Arm Rest Assy w/ Arm Pad		\$40.00**
<input type="checkbox"/>	104400-78700	(R) Arm Rest Assy w/ Arm Pad		\$40.00**

TIRES

<input type="checkbox"/>	TSR20050	200X50 Solid Tire Rear	x2	\$88.00**
<input type="checkbox"/>	105101-78700	Solid Front Wheel Assy-7x1 1/4	x2	\$180.00**

ANTI-TIP WHEELS

<input type="checkbox"/>	209928-89300	Anti-Tip Wheel Assy	x2	\$40.00**
--------------------------	--------------	---------------------	----	-----------

CHARGER

<input type="checkbox"/>	109201-88801-A2	2 amp Charger Off-Board		\$130.00**
--------------------------	-----------------	-------------------------	--	------------

ACCESSORIES

<input type="checkbox"/>	209916-78700	Wire Mesh Front Basket		\$40.00**
<input type="checkbox"/>	OXY-03	Mesh Oxygen Tank Holder		\$60.00
<input type="checkbox"/>	WIG-0168	Wig-wag extensions (pair)		\$20.00
<input type="checkbox"/>	300909-03	Safety Belt		\$34.00

SCOOTIE SPECIFICATIONS

OVERALL DIMENSIONS

Base Length.....	37"-40"
Base Width.....	17"
Height.....	34"
(Ground to Seat Top)	
Ground Clearance.....	2"

WEIGHTS

User Weight Capacity.....	250 Lbs
Total Unit Weight.....	82 Lbs
Weight of Heaviest Part.....	25 Lbs
Battery Module Weight (Pack).....	22 Lbs

TIRES

Number, Size of Front Tire.....	2, 7"
Number, Size of Rear Tire.....	2, 8"

BATTERIES

BAT-12 Module Capacity.....	12v12Ah x2
-----------------------------	------------

CHARGER

Charger-On Board.....	2 amp
-----------------------	-------

MOTOR POWER OUTPUT

Horsepower.....	.75
-----------------	-----

DRIVING SPECS

Turning Radius.....	38"
Range.....	10 miles
Maximum Speed.....	3.5 mph
Maximum Climbing Angle.....	14% (8°)

TOTAL \$

****STANDARD WITH UNIT**

All specifications are subject to change without notice.

Shoprider Mobility Products, Inc® reserves the rights of any changes to the unit:

- 1) Include the anti-tip wheel or rear castor.
- 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown below:
 - a) The weight of user.
 - b) Ground surface and terrain.
 - c) Battery capacity and conditions.
 - d) Type of charger.
 - e) Ambient Temperature.
 - f) The way and/or method of driving.
 - g) If the battery and mechanical moving parts are fully broken in.