

# PRICE LIST & ORDER FORM

# SHOPRIDER®

For easy ordering and information, please contact Customer Service **Toll Free at 1.800.743.0772 FAX: 1.877.797.7081** Effective July 2009



**TE-787 NA SCOOTIE**

NAME: \_\_\_\_\_  
 PHONE NO: \_\_\_\_\_  
 ACCT NO: \_\_\_\_\_  
 PO #: \_\_\_\_\_ DATE: \_\_\_\_\_  
 ADDRESS: \_\_\_\_\_  
 CITY: \_\_\_\_\_ STATE: \_\_\_\_\_ ZIP: \_\_\_\_\_  
 SHIP TO: \_\_\_\_\_  
 PHONE NO: \_\_\_\_\_  
 ADDRESS: \_\_\_\_\_  
 CITY: \_\_\_\_\_ STATE: \_\_\_\_\_ ZIP: \_\_\_\_\_

## SPECIFICATIONS

### OVERALL DIMENSIONS

Base Length.....37"  
 Base Width..... 17"  
 Height.....30.5"  
 Ground to Seat Top  
 Ground Clearance.....2"

### WEIGHTS

User Weight Capacity.....250 Lb  
 Total Unit Weight.....82 Lbs  
 Weight of Heaviest Part.....25 Lbs  
 Battery Weight Module (pack) .....22 Lbs

### TIRES

Number, Size of Front Tire.....2, 7"  
 Number, Size of Rear Tire.....2, 8"

### BATTERIES

BAT -12 Module Capacity.....12v12Ah x2

### CHARGER

Charger-On Board.....2 amp

### MOTOR POWER OUTPUT

Horsepower.....0.75 hp

### DRIVING SPECS

Turning Radius .....38"  
 Range .....10 miles  
 Maximum Speed .....3.5 mph  
 Maximum Climbing Angle .....14% / 8°

All specifications are subject to change without notice. Shoprider Mobility Products, Inc.® reserves the rights of any changes to the unit:

- 1) Include the anti-tip wheel or rear castor.
- 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown below:
  - a) The weight of user.
  - b) Ground surface and terrain.
  - c) Battery capacity and conditions.
  - d) Type of charger.
  - e) Ambient Temperature.
  - f) The way and/or method of driving.
  - g) If the battery and mechanical moving parts are fully broken in.

### BASE SCOOTER

<input type="checkbox"/> TE-787NA	SCOOTIE RED	K0800	\$1,595.00
<input type="checkbox"/> TE-787NA	SCOOTIE BLUE	K0800	\$1,595.00

### BATTERY

<input type="checkbox"/> BAT-12	12V 12AMP	x2	\$75.00**ea
---------------------------------	-----------	----	-------------

### SEATS (W"xD"xH")

<input type="checkbox"/> SF151413G	Grey Vinyl Fish-On Seat		\$100.00**
------------------------------------	-------------------------	--	------------

### ARM REST ASSY

<input type="checkbox"/> 104300-78700	(L) Arm Rest Assy w/ Arm Pad		\$40.00**
<input type="checkbox"/> 104400-78700	(R) Arm Rest Assy w/ Arm Pad		\$40.00**

### TIRES

<input type="checkbox"/> TSR20050	200X50 Solid Tire Rear	x2	\$44.00**ea
<input type="checkbox"/> 105101-78700	Solid Front Wheel Assy-7x1¼	x2	\$90.00**ea

### ANTI-TIP WHEELS

<input type="checkbox"/> 209928-89300	Anti-Tip Wheel Assy	x2	\$40.00**
---------------------------------------	---------------------	----	-----------

### CHARGER

<input type="checkbox"/> 109201-88801-A2	2 amp Charger Off-Board		\$130.00**
--	-------------------------	--	------------

### ACCESSORIES

<input type="checkbox"/> 209916-78700	Wire Mesh Front Basket		\$40.00**
<input type="checkbox"/> OXY-03	Mesh Oxygen Tank Holder		\$60.00
<input type="checkbox"/> WIG-0168	Wig-wag extensions (pair)		\$20.00
<input type="checkbox"/> 300909-03	Safety Belt		\$34.00

\*\*STANDARD WITH UNIT

TOTAL \$ \_\_\_\_\_