



PRICE LIST & ORDER FORM Effective November 2006

For easy ordering and information contact Customer Service Toll Free at:
800-743-0772 Fax: 877-797-7081

Name _____	Phone No _____
Acct No _____ PO # _____	Date _____
Address _____	City/State _____ Zip _____
Ship To _____	Phone No _____
Address _____	City/State _____ Zip _____
Comments _____	PO # _____

Stock For your convenience, order by stock number (e.g. BAT-U1)

BASE CHAIR

888WHD-BGRD SHOPRIDER HD /Burgundy K0827 \$9,650.00

BATTERY

BAT-75 12V 75AMP x2 \$250.00ea.

SEATS (W"xD"xH")

SV222120RG 22x21x20 EXTRA WIDE, ADJUSTABLE SEAT \$420.00 **

ARM REST ASSY

1343-00205-00 (R) Adj Ht / Wd Swing-up Arm Rest Assy \$200.00 **

1343-00105-00 (L) Adj Ht / Wd Swing-up Arm Rest Assy \$200.00 **

204304-001403 ARM PAD \$84.00 **

TIRES

TSR26085 10" Flat-Free Rear Tires x2 \$50.00ea. **

CASTORS

TSF280250 9"Front Flat Free Castor Tire x2 \$50.00ea. **

1315-00107-00 Castor Front Fork x2 \$128.00ea. **

ANTI-TIP WHEELS

800405-001-00 R) Anti-tip \$72.00 **

800405-001-01 (L) Anti-tip \$72.00 **

FOOT REST

800907-007-06 Ex-Large Flip Up Foot Plate Assy \$560.00 **

ELRHD-B(K) ELR Kit (Elevating Leg Rests-Black) \$340.00

800907-010-02 Hangers, Receivers, Hardware set+ \$188.00

1329-00403-00 Swing Away Foot Rest - Left \$80.00

1329-00404-00 Swing Away Foot Rest - Right \$80.00

CONTROLLER

1361-00114-00 VR-2 Remote Joystick Controller \$72.00 **

1361-00104-00 VR-2 Main Control Box \$776.00 **

CHARGER

109201-88500 8 amp Charger Off-Board \$300.00 **

ACCESSORIES

CCH-003(K) Cane Holder \$150.00

CH-02(K) Oxygen Holder (Metal) \$150.00

WH-0185(K) Walker Holder \$150.00

OXY-03(K) Mesh Oxygen Tank Holder \$60.00

300909-03-B1 Safety Belt 66" Long \$40.00 **

CUPHOLDER(K) CUP HOLDER \$30.00

RB-0125(K) Rear Basket Large (21"W15"L9"D) \$170.00

TOTAL \$

SHOPRIDER HD SPECIFICATIONS

OVERALL DIMENSIONS

Base Length.....47"

Base Width.....29"

Height.....50"

(Ground to Seat Top)

Ground Clearance.....3"

WEIGHTS

User Weight Capacity.....500 Lbs

Total Unit Weight.....304 Lbs

Weight of Heaviest Part.....n/a

Battery Module Weight (75amp)....57lbsx2

TIRES/CASTORS

Front Castor Assy Flat Free.....9"

Rear Flat-Free Tires.....10"

BATTERIES

Battery Module Capacity.....12v75Ah x2

CHARGER

Charger-Off Board.....8 amp

MOTOR POWER OUTPUT

Horsepower.....1hp x 2

DRIVING SPECS

Turning Radius.....37.5"

Range.....up to 31 miles

Maximum Speed.....4.375mph

Maximum Climbing Angle.....10% / 6°

All specifications are subject to change without notice.

Shoprider Mobility Products, Inc® reserves the rights of any changes to the unit:

1) Include the anti-tip wheel or rear castor.

2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.

3) The actual driving range varies with the factors shown below:

a) The weight of user.

b) Ground surface and terrain.

c) Battery capacity and conditions.

d) Type of charger.

e) Ambient Temperature.

f) The way and/or method of driving.

g) If the battery and mechanical moving parts are fully broken in.