

# SHOPPRIDER®



## UL8W SMARTIE

### PRICE LIST & ORDER FORM Effective July 2009

For easy ordering and information contact Customer Service Toll Free at:

**800-743-0772 Fax: 877-797-7081**

Name _____	Phone No _____
Acct No _____	PO # _____
Address _____	City/State _____ Zip _____
Ship To _____	Phone No _____
Address _____	City/State _____ Zip _____

stock #	item description		msrp
<b>BASE SCOOTER</b>			
<input type="checkbox"/> UL8W	Smartie Powerchair	K0821	\$5,495.00
<b>BATTERY</b>			
<input type="checkbox"/> PH012010D01	Red or Blue Envirofriendly Battery Pk	x1	\$740.00**
<i>***Envirofriendly Battery must be purchased with Smartie</i>			
<b>SEATS (W"xD"xH")</b>			
<input type="checkbox"/> SV1817½19G	18" Grey Vinyl Deluxe Captain Seat		\$400.00**
<b>ARM REST ASSY</b>			
<input type="checkbox"/> 104300-00100-B2	(L) Arm Rest Assy w/ Arm Pad		\$160.00**
<input type="checkbox"/> 104400-00100-B2	(R) Arm Rest Assy w/ Arm Pad		\$160.00**
<b>TIRES</b>			
<input type="checkbox"/> 105201-69101	200X70 Solid Rear Wheel Assy	x2	\$90.00**ea
<input type="checkbox"/> 1351-00601-00	6" Solid Tire Front	x2	\$140.00**ea
<b>ANTI-TIP WHEELS</b>			
<input type="checkbox"/> 209928-89300	Anti-Tip Wheel Assy	x2	\$20.00**ea
<b>CHARGER</b>			
<input type="checkbox"/> 109201-88805	3 pin; 2 Amp charger; UL8W; Li		\$100.00**
<b>CONTROLLER</b>			
<input type="checkbox"/> 1361-00300-01	Joystick Controller-UL8WPB		\$998.00**
<b>FOOT REST</b>			
<input type="checkbox"/> 2329-00100-10	Full Function Foot Plate		\$200.00**
<input type="checkbox"/> ELRUL-R(K)	Elevating Leg Rest Kit in Red		\$335.00
<input type="checkbox"/> ELRUL-B(K)	Elevating Leg Rest Kit in Black		\$335.00
<b>ACCESSORIES</b>			
<input type="checkbox"/> OXY-03	Mesh Oxygen Tank Holder		\$60.00
<input type="checkbox"/> CCH-003-GK(K)	Cane Holder Kit		\$150.00
<input type="checkbox"/> WH-0185-GK(K)	Walker Holder Kit		\$150.00

SPECIFICATIONS	
<b>OVERALL DIMENSIONS</b>	
Base Length.....	34.3"
w/o footplate.....	25.3"
Base Width.....	21.5"
Height.....	30.5"
(Ground to Seat Top)	
Ground Clearance.....	2"
<b>WEIGHTS</b>	
User Weight Capacity.....	250 Lbs
Front Chassis w. batteries.....	27.1Lbs.
Rear Chassis w/ drive train.....	29 Lbs.
Total Weight(w/ Bat's & Seat).....	93 Lbs
Battery pack weight.....	5.9 Lbs.
Captain seat w/controller assy....	37.1Lbs
Battery Module Weight (Pack).....	5.9 Lbs
<b>TIRES</b>	
Number, Size of Front Tire.....	2, 6"
Number, Size of Rear Tire.....	2, 8"
<b>BATTERIES</b>	
Envirofriendly Battery Pack.....	5.9 Lbs
<b>CHARGER</b>	
Charger-Off Board.....	2 amp
<b>MOTOR POWER OUTPUT</b>	
Horsepower.....	.45hp X2
<b>DRIVING SPECS</b>	
Turning Radius.....	15.5"
Range.....	10 miles
Maximum Speed.....	3.75 mph
Maximum Climbing Angle.....	10% (6°)

**TOTAL \$**

**\*\*STANDARD equipment**

All specifications are subject to change without notice.  
Shopprider Mobility Products, Inc® reserves the rights of any changes to the unit:

- 1) Include the anti-tip wheel or rear castor.
- 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown below:
  - a) The weight of user.
  - b) Ground surface and terrain.
  - c) Battery capacity and conditions.
  - d) Type of charger.
  - e) Ambient Temperature.
  - f) The way and/or method of driving.
  - g) If the battery and mechanical moving parts are fully broken in.