



**GK9**  
**START 8**

**PRICE LIST & ORDER FORM** Effective May 2007

For easy ordering and information contact Customer Service Toll Free at:  
**800-743-0772 Fax: 877-797-7081**

Name _____	Phone No _____
Acct No _____	PO # _____
Address _____	Date _____
City/State _____	Zip _____
Ship To _____	Phone No _____
Address _____	City/State _____
	Zip _____

Stock	For your convenience, order by stock number (e.g. BAT-1/2)		
<b>BASE SCOOTER</b>			
<input type="checkbox"/>	GK9	START 8 Burgundy	\$2,395.00
<b>BATTERY</b>			
<input type="checkbox"/>	please call	Battery Pack w/ 2 x 12v-21amp SLA Batteries	xxxxx
<input type="checkbox"/>	BAT-1/2	12V - 21amp SLA battery x2	\$90.00 ea**
		*Complimentary Battery has limited warranty	
<b>SEATS (W"xD"xH")</b>			
<input type="checkbox"/>	please call	18"Wx16"D Stadium Style Fish-On Assy	xxxxx
<input type="checkbox"/>	please call	Seat/Back Cushion Set-Vinyl Pads	xxxxx
<input type="checkbox"/>	SP161618RG	Upgrade 16"w Seat Pan/Capt Back Recline*	\$150.00
		*Also avail in 18"w. Seat cushion not sold.	
<b>ARM REST ASSY</b>			
<input type="checkbox"/>	104300-00100-B2	(L) Flip-up Arm Rest Assy w/ Arm Pad	\$160.00**
<input type="checkbox"/>	104400-00100-B2	(R) Flip-up Arm Rest Assy w/ Arm Pad	\$160.00**
<input type="checkbox"/>	104310-00102	Arm pad w/ height-adjustment tube	\$100.00**
<b>TIRES</b>			
<input type="checkbox"/>	please call	Solid Tire Rear Wheel Assy x2	xxxxx
<input type="checkbox"/>	please call	Solid Tire Front Wheel Assy x2	xxxxx
<b>ANTI-TIP WHEELS</b>			
<input type="checkbox"/>	please call	Anti-Tip Wheel Assy x2	xxxxx
<b>CHARGER</b>			
<input type="checkbox"/>	109201-88306	3 amp Charger Off-Board	\$180.00**
<b>ACCESSORIES</b>			
<input type="checkbox"/>	209916-88100	Wire Mesh Front Basket	\$50.00**
<input type="checkbox"/>	OXY-03	Mesh Oxygen Tank Holder	\$60.00
<input type="checkbox"/>	CH-02-GK(K)	Metal Oxygen Holder Kit	\$150.00
<input type="checkbox"/>	CCH-003-GK(K)	Cane Holder Kit	\$150.00
<input type="checkbox"/>	RB-0125-GK(K)	Rear basket Kit	\$150.00
<input type="checkbox"/>	WH-0185-GK(K)	Walker Holder Kit	\$150.00
<input type="checkbox"/>	WIG-0168	Wig-wag extensions (pair)	\$20.00
<input type="checkbox"/>	300909-03	Safety Belt	\$34.00

START 8 SPECIFICATIONS	
<b>OVERALL DIMENSIONS</b>	
Base Length.....	42"
Base Width.....	22"
Height.....	19-21"
(Ground to Seat Top)	
Min. Obstacle Height.....	1.5"
<b>WEIGHTS</b>	
User Weight Capacity.....	300 Lbs
Total Unit Weight.....	123 Lbs
Weight of Heaviest Part.....	36 Lbs
Battery Module Weight (Pack).....	32 Lbs
<b>TIRES</b>	
Number, Size of Front Tire.....	2, 9"
Number, Size of Rear Tire.....	2, 9"
<b>BATTERIES</b>	
Batt Pk Module Capacity.....	12v21Ah x2
<b>CHARGER</b>	
Charger-Off Board.....	3 amp
<b>MOTOR POWER OUTPUT</b>	
Horsepower.....	.75
<b>DRIVING SPECS</b>	
Turning Radius.....	45.3"
Range.....	up to 12.5 mis
Maximum Speed.....	3.8 mph
Maximum Climbing Angle.....	10% (8°)
at 250 lbs user weight	

**TOTAL \$**

**\*\*STANDARD WITH UNIT**

- All specifications are subject to change without notice.  
Shoprider Mobility Products, Inc® reserves the rights of any changes to the unit:
- 1) Include the anti-tip wheel or rear castor.
  - 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
  - 3) The actual driving range varies with the factors shown below:
    - a) The weight of user.
    - b) Ground surface and terrain.
    - c) Battery capacity and conditions.
    - d) Type of charger.
    - e) Ambient Temperature.
    - f) The way and/or method of driving.
    - g) If the battery and mechanical moving parts are fully broken in.