



**GK9
START 8**

PRICE LIST & ORDER FORM Effective May 2007

For easy ordering and information contact Customer Service Toll Free at:
800-743-0772 Fax: 877-797-7081

Name _____	Phone No _____
Acct No _____	PO # _____
Address _____	City/State _____ Zip _____
Ship To _____	Phone No _____
Address _____	City/State _____ Zip _____

Stock	For your convenience, order by stock number (e.g. BAT-1/2)		
BASE SCOOTER			
<input type="checkbox"/>	GK9	START 8 Burgundy	\$2,395.00
BATTERY			
<input type="checkbox"/>	800501-707-00	Battery Pack w/ 2 x 12v-21amp SLA Batteries	\$312.00
<input type="checkbox"/>	BAT-21	12V - 21amp SLA battery x2	\$110.00 ea**
*Complimentary Battery has limited warranty			
SEATS (W"xD"xH")			
<input type="checkbox"/>	please call	18"Wx16"D Stadium Style Fish-On Assy	xxxxx
<input type="checkbox"/>	please call	Seat/Back Cushion Set-Vinyl Pads	xxxxx
<input type="checkbox"/>	SP161618RG	Upgrade 16"w Seat Pan/Capt Back Recline*	\$150.00
*Also avail in 18"w. Seat cushion not sold.			
ARM REST ASSY			
<input type="checkbox"/>	104300-00100-B2	(L) Flip-up Arm Rest Assy w/ Arm Pad	\$160.00**
<input type="checkbox"/>	104400-00100-B2	(R) Flip-up Arm Rest Assy w/ Arm Pad	\$160.00**
<input type="checkbox"/>	104310-00102	Arm pad w/ height-adjustment tube	\$100.00**
TIRES			
<input type="checkbox"/>	please call	Solid Tire Rear Wheel Assy x2	xxxxx
<input type="checkbox"/>	please call	Solid Tire Front Wheel Assy x2	xxxxx
ANTI-TIP WHEELS			
<input type="checkbox"/>	please call	Anti-Tip Wheel Assy x2	xxxxx
CHARGER			
<input type="checkbox"/>	109201-88306	3 amp Charger Off-Board	\$230.00**
ACCESSORIES			
<input type="checkbox"/>	209916-88100	Wire Mesh Front Basket	\$50.00**
<input type="checkbox"/>	OXY-03	Mesh Oxygen Tank Holder	\$60.00
<input type="checkbox"/>	CH-02-GK(K)	Metal Oxygen Holder Kit	\$150.00
<input type="checkbox"/>	CCH-003-GK(K)	Cane Holder Kit	\$150.00
<input type="checkbox"/>	RB-0125-GK(K)	Rear basket Kit	\$150.00
<input type="checkbox"/>	WH-0185-GK(K)	Walker Holder Kit	\$150.00
<input type="checkbox"/>	WIG-0168	Wig-wag extensions (pair)	\$20.00
<input type="checkbox"/>	300909-03	Safety Belt	\$34.00

START 8 SPECIFICATIONS	
OVERALL DIMENSIONS	
Base Length.....	42"
Base Width.....	22"
Height.....	19-21"
(Ground to Seat Top)	
Min. Obstacle Height.....	1.5"
WEIGHTS	
User Weight Capacity.....	300 Lbs
Total Unit Weight.....	123 Lbs
Weight of Heaviest Part.....	36 Lbs
Battery Module Weight (Pack).....	32 Lbs
TIRES	
Number, Size of Front Tire.....	2, 9"
Number, Size of Rear Tire.....	2, 9"
BATTERIES	
Batt Pk Module Capacity.....	12v21Ah x2
CHARGER	
Charger-Off Board.....	3 amp
MOTOR POWER OUTPUT	
Horsepower.....	.75
DRIVING SPECS	
Turning Radius.....	45.3"
Range.....	up to 12.5 mis
Maximum Speed.....	3.8 mph
Maximum Climbing Angle.....	10% (8°)
at 250 lbs user weight	

TOTAL \$

****STANDARD WITH UNIT**

- All specifications are subject to change without notice.
Shoprider Mobility Products, Inc® reserves the rights of any changes to the unit:
- 1) Include the anti-tip wheel or rear castor.
 - 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
 - 3) The actual driving range varies with the factors shown below:
 - a) The weight of user.
 - b) Ground surface and terrain.
 - c) Battery capacity and conditions.
 - d) Type of charger.
 - e) Ambient Temperature.
 - f) The way and/or method of driving.
 - g) If the battery and mechanical moving parts are fully broken in.