

PRICE LIST & ORDER FORM

SHOPRIDER®

For easy ordering and information, please contact Customer Service **Toll Free at 1.800.743.0772 FAX: 1.877.797.7081** Effective July 2010



888WA STREAMER SPORT

NAME: _____
 PHONE NO: _____
 ACCT NO: _____
 PO #: _____ DATE: _____
 ADDRESS: _____
 CITY: _____ STATE: _____ ZIP: _____
 SHIP TO: _____
 PHONE NO: _____
 ADDRESS: _____
 CITY: _____ STATE: _____ ZIP: _____

SPECIFICATIONS

OVERALL DIMENSIONS

Base Length.....41"
 Base Width.....23"
 Height.....40"
 Ground to Seat Top
 Ground Clearance.....3"

WEIGHTS

User Weight Capacity.....300 Lbs
 Total Unit Weight.....164 Lbs
 Battery Weight Module (U1).....22.5 Lbs

TIRES/ CASTORS

Front castor Assy.....6"
 Rear Flat Free.....10"

BATTERIES

BAT-U1 Module Capacity.....12v32Ah x2

CHARGER

Charger-Off Board.....5 amp

MOTOR POWER OUTPUT

Horsepower.....0.75 x2

DRIVING SPECS

Turning Radius17.7"
 Rangeup to 10 miles
 Maximum Speed5 mph
 Maximum Climbing Angle14% / 8°

All specifications are subject to change without notice. Shoprider Mobility Products, Inc® reserves the rights of any changes to the unit:

- 1) Include the anti-tip wheel or rear castor.
- 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown below:
 - a) The weight of user.
 - b) Ground surface and terrain.
 - c) Battery capacity and conditions.
 - d) Type of charger.
 - e) Ambient Temperature.
 - f) The way and/or method of driving.
 - g) If the battery and mechanical moving parts are fully broken in.

BASE CHAIR

888WA-RED Streamer WA (prejudgment) K0823 \$5,295.00
 with seat pan upgrade K0822

BATTERY

BAT-U1 12V 32AMP x2 \$140.00ea**

SEATS (W"xD"xH")

SV191821B 19X18X21 Captain Seat Black \$280.00**
 SV191820RG 19x18x20 Captain Seat Grey, w/ Reclining \$400.00
 SV201820RB 20x18x20 Black vinyl leatherette w. pillow-top \$370.00
 SV2121211/2RB 21x21x21 1/2 Captain seat; black vinyl; 21" \$400.00
 SV2121211/2RG 21x21x21 1/2 Captain seat; grey vinyl; 21" \$400.00
 SP181821RG 18x18x21 Pan seat; 2-tone grey (18x18) \$400.00
 SP202023RG 20X20X23 Pan seat; 2-tone grey (20x20) \$400.00
 SP161618RG 16x16x18 Pan seat; 2-tone grey reclining \$400.00
 Seat Cushion NOT INCLUDED when ordering Seat Pan

ARM REST ASSY

104300-00106-A3 (L) Flipback High Adjustable Armrest \$140.00**
 104400-00106-A3 (R) Flipback High Adjustable Armrest \$140.00**

FRONT CASTORS

1351-00601-00 6" Castor Assy w/fork x2 \$140.00ea**
 205201-00600-A2 6" Castow w/o fork x2 \$88.00 ea

DRIVE TIRES

TSR26085 10" Flat-Free Tires (260x85) x2 \$50.00ea**
 TAR26085 10" Pnuematic Tires (260x85) x2 \$32.00ea
 205102-88100 10" Tire Tubes (260x85) x2 \$20.00ea

FOOT REST

2329-00100-11 Flip Up Foot Plate \$200.00**
 ELR-STREAMER ELR Kit complete Streamer- 888WA \$180.00

CONTROLLER

1361-00114-00 VR2 Joystick Controller \$360.00**
 1361-00104-00 VR2 Main Control Box \$776.00**

CHARGER

109201-88124 5 amp Charger Off-Board \$230.00**

ACCESSORIES

BKS100L(K) Left Stump Support \$192.00
 BKS100R(K) Right Stump Support \$192.00
 1343-00709-00 Right Swing Away Mount Bracket for VSI \$120.00
 1343-00709-01 Left Swing Away Mount Bracket for VSI \$120.00
 CCH-003(K) Cane Holder \$150.00
 O210(K) Oxygen Tank Holder (Metal) \$150.00
 WH-0185(K) Walker Holder \$150.00
 OXY-03 Mesh Oxygen Tank Holder \$60.00
 300909-03 Safety Belt \$34.00**
 P75755 VSI Joystick Knob \$12.00**

**STANDARD WITH UNIT

TOTAL \$ _____