



**777-3/4
SUNRUNNER S**

PRICE LIST & ORDER FORM Effective April 2006

For easy ordering and information contact Customer Service Toll Free at:
800-743-0772 Fax: 877-797-7081

Name _____	Phone No _____
Acct No _____	PO # _____
Address _____	City/State _____ Zip _____
Ship To _____	Phone No _____
Address _____	City/State _____ Zip _____

Stock For your convenience, order by stock number (e.g. BAT-U1)

BASE SCOOTER

<input type="checkbox"/>	777-3/4	SUNRUNNER S	BURGUNDY	E1230	\$1,995.00
<input type="checkbox"/>	777-3/4	SUNRUNNER S	BLUE	E1230	\$1,995.00

BATTERY

<input type="checkbox"/>	BAT-12	12V-12amp SLA Battery	x2	\$70.00 ea**
--------------------------	--------	-----------------------	----	--------------

COMPLIMENTARY BATTERIES ARE NOT WARRANTED

SEATS (W"xD"xH")

<input type="checkbox"/>	SF181614G	18"w Fish-On Seat w/ Seat Pads-Vinyl	\$218.00**
<input type="checkbox"/>	SV171714RB	Upgrade 17"x17"x14" Recline Blk Pillow-Top	\$150.00
<input type="checkbox"/>	SP161618RG	Upgrade 16" Seat Pan/Capt Back Recline *Seat Cushion not included	\$150.00

ARM REST ASSY

<input type="checkbox"/>	104400-88800	(L) Arm Rest Assy w/ Arm Pad	\$160.00**
<input type="checkbox"/>	104300-88800	(R) Arm Rest Assy w/ Arm Pad	\$160.00**

TIRES

<input type="checkbox"/>	TSF20050	200x50 Solid Tire-Front	x2	\$50.00 ea**
<input type="checkbox"/>	TSR280250	2.80x2.50 Solid Tire-Rear	x2	\$56.00 ea**

ANTI-TIP WHEELS

<input type="checkbox"/>	209928-89300	Anti-Tip Wheel Assy	x2	\$20.00**
--------------------------	--------------	---------------------	----	-----------

CHARGER

<input type="checkbox"/>	109201-88801-A2	2 amp Charger Off-Board	\$100.00**
--------------------------	-----------------	-------------------------	------------

ACCESSORIES

<input type="checkbox"/>	209916-78700	Wire Mesh Front Basket	\$40.00**
<input type="checkbox"/>	CCH-003(K)	Cane Holder	\$150.00
<input type="checkbox"/>	CH-02(K)	Oxygen Tank Holder (Metal)	\$150.00
<input type="checkbox"/>	WH-0185(K)	Walker Holder	\$150.00
<input type="checkbox"/>	OXY-03	Mesh Oxygen Tank Holder	\$60.00
<input type="checkbox"/>	WIG-0168	Wig-wag extensions (pair)	\$20.00
<input type="checkbox"/>	RB-0125(K)	Rear Wire Mesh Basket	\$170.00
<input type="checkbox"/>	300909-03	Safety Belt	\$34.00

TOTAL \$ _____

SUNRUNNER S SPECIFICATIONS	
OVERALL DIMENSIONS	
Base Length.....	43"-46"
Base Width.....	22"
Height.....	35"
(Ground to Seat Top)	
Ground Clearance.....	3"
WEIGHTS	
User Weight Capacity.....	250 Lbs
Total Unit Weight.....	139 Lbs
Weight of Heaviest Part.....	40 Lbs
Battery Module Weight (pack).....	22 Lbs
TIRES	
Number, Size of Front Tire.....	2, 8"
Number, Size of Rear Tire.....	2, 9"
BATTERIES	
BAT-12 Module Capacity.....	12v12Ah x2
CHARGER	
Charger-Off Board.....	2 amp
MOTOR POWER OUTPUT	
Horsepower.....	.75
DRIVING SPECS	
Turning Radius.....	37"
Range.....	10 miles
Maximum Speed.....	5 mph
Maximum Climbing Angle.....	14% (8°)

****STANDARD WITH UNIT**

All specifications are subject to change without notice.

Shoprider Mobility Products, Inc® reserves the rights of any changes to the unit:

- 1) Include the anti-tip wheel or rear castor.
- 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown below:
 - a) The weight of user.
 - b) Ground surface and terrain.
 - c) Battery capacity and conditions.
 - d) Type of charger.
 - e) Ambient Temperature.
 - f) The way and/or method of driving.
 - g) If the battery and mechanical moving parts are fully broken in.