



PRICE LIST & ORDER FORM



For easy ordering and information, please contact Customer Service **Toll Free at 1.800.743.0772 FAX: 1.877.797.7081** Effective July 2009

 <p>888B-3 SUNRUNNER 3</p>	 <p>888B-4 SUNRUNNER 4</p>	NAME: _____
		PHONE NO: _____
		ACCT NO: _____
		PO #: _____ DATE: _____
		ADDRESS: _____
		CITY: _____ STATE: _____ ZIP: _____
		SHIP TO: _____
		PHONE NO: _____
ADDRESS: _____		
CITY: _____ STATE: _____ ZIP: _____		

SPECIFICATIONS	
SUNRUNNER-3 (888B-3)	
Base LxWxH*.....	44"x23"x42"
*Ground to Seat Top	
Ground Clearance.....	3"
User Weight Capacity.....	300 lbs
Total Unit Weight.....	165 lbs
Battery Module Weight (U1).....	25 lbs x2
Number, Size Tire (front/rear).....	1, 9"/ 2, 10"
BAT-U1 Module Capacity.....	12V35Ahx2
Charger-Off Board.....	3 amp
Turning Radius.....	37"
Range.....	up to 25 miles
Maximum Speed.....	5 mph
Maximum Climbing Angle.....	14% /8°
SUNRUNNER-4 (888B-4)	
Base LxWxH*.....	48"x23"x43"
*Ground to Seat Top	
Ground Clearance.....	3"
User Weight Capacity.....	300 lbs
Total Unit Weight.....	180 lbs
Battery Module Weight (U1).....	25 lbs x2
Number, Size Tire (front/rear).....	2, 10"/ 2, 10"
BAT-U1 Module Capacity.....	12V35Ahx2
Charger-Off Board.....	3 amp
Turning Radius.....	50"
Range.....	up to 25 miles
Maximum Speed.....	5 mph
Maximum Climbing Angle.....	14% /8°

All specifications are subject to change without notice. Shoprider Mobility Products, Inc.® reserves the rights of any changes to the unit:
 1) Include the anti-tip wheel or rear castor.
 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
 3) The actual driving range varies with the factors shown below:
 a) The weight of user.
 b) Ground surface and terrain.
 c) Battery capacity and conditions.
 d) Type of charger.
 e) Ambient Temperature.
 f) The way and/or method of driving.
 g) If the battery and mechanical moving parts are fully broken in.

BASE SCOOTER		MSRP
<input type="checkbox"/> 888B-3	SUNRUNNER-3 Burgundy/Blue	K0806 \$2,195.00
<input type="checkbox"/> 888B-4	SUNRUNNER-4 Burgundy/Blue	K0899 \$2,395.00
BATTERY		
<input type="checkbox"/> BAT-U1	12V 35AMP	x2 \$110.00**
SEATS (W"xD"xH")		
<input type="checkbox"/> SV171/21815G	17½" Grey Vinyl Deluxe Captain Seat	\$408.00**
<input type="checkbox"/> SV171/21714RB	Upgrade 17½"x17" Pillow-Top Recline Seat Vinyl	\$150.00
<input type="checkbox"/> SV201820RG	Upgrade 20"x18" Recline Capt Grey Vinyl	\$150.00
<input type="checkbox"/> SP181821RG	Upgrade Pan seat w/ Capt back recline (18x18)	\$150.00
ARM REST ASSY		
<input type="checkbox"/> 104300-88800	(L) Arm Rest Assy w/ Arm Pad	\$160.00**
<input type="checkbox"/> 104400-88800	(R) Arm Rest Assy w/ Arm Pad	\$160.00**
TIRES		
<input type="checkbox"/> TAR26085 Rear	10" Pnuematic Tires-all Models	x2 \$32.00**
<input type="checkbox"/> 205102-88100 Rear	10" Tire Tubes-all Models	x2 \$20.00**
<input type="checkbox"/> TAF26085 Front	10" Pnuematic Tires-888B-4	x2 \$32.00**
<input type="checkbox"/> 205102-88100 Front	10" Tire Tubes-888B-4	x2 \$20.00**
<input type="checkbox"/> TAF280250 Front	9" Pnuematic Tires-888B-3	x1 \$36.00**
<input type="checkbox"/> 205202-77100 Front	9" Tire Tubes-888B-3	x1 \$20.00**
<input type="checkbox"/> STU8-9(K) Front & Rear	888B-3 Solid Tire Upgrade Kit	\$150.00
<input type="checkbox"/> STU8A(K) Front & Rear	888B-4 Solid Tire Upgrade Kit	\$200.00
CHARGER		
<input type="checkbox"/> 109201-88124	5 amp Charger Off-Board	\$230.00**
<input type="checkbox"/> 109201-88306	3 amp Charger Off-Board	\$230.00
ACCESSORIES		
<input type="checkbox"/> 109901-88101	Front Wicker Basket	\$70.00**
<input type="checkbox"/> CCH-003(K)	Cane Holder	\$150.00
<input type="checkbox"/> 0210(K)	Oxygen Tank Holder (Metal)	\$150.00
<input type="checkbox"/> WH-0185(K)	Walker Holder	\$150.00
<input type="checkbox"/> WIG-0168	Wig-wag extensions (pair)	\$20.00
<input type="checkbox"/> RB-0125(K)	Rear Wire Mesh Basket	\$170.00
<input type="checkbox"/> 300909-03	Safety Belt	\$34.00

**STANDARD WITH UNIT

TOTAL \$ _____