



**888B-3**  
Sunrunner-3



**888B-4**  
Sunrunner-4

**PRICE LIST & ORDER FORM** Effective Novembe 2006

For easy ordering and information contact  
Customer Service Toll Free at:

**800-743-0772 Fax: 877-797-7081**

|               |                            |
|---------------|----------------------------|
| Name _____    | Phone No _____             |
| Acct No _____ | PO # _____                 |
| Date _____    |                            |
| Address _____ | City/State _____ Zip _____ |
| Ship To _____ | Phone No _____             |
| Address _____ | City/State _____ Zip _____ |

**Stock** For your convenience, order by stock number (e.g. BAT-U1)

| BASE SCOOTER             |        |             |               |       |            |
|--------------------------|--------|-------------|---------------|-------|------------|
| <input type="checkbox"/> | 888B-3 | SUNRUNNER-3 | Burgundy/Blue | K0806 | \$2,195.00 |
| <input type="checkbox"/> | 888B-4 | SUNRUNNER-4 | Burgundy/Blue | K0899 | \$2,395.00 |

| BATTERY                  |        |           |  |    |            |
|--------------------------|--------|-----------|--|----|------------|
| <input type="checkbox"/> | BAT-U1 | 12V 35AMP |  | x2 | \$110.00** |

COMPLIMENTARY BATTERIES ARE NOT WARRANTED

| SEATS (W"xD"xH")         |               |  |  |  |            |
|--------------------------|---------------|--|--|--|------------|
| <input type="checkbox"/> | SV1717131/2G  | 17½" Deluxe Captain Seat Vinyl-888B-3 & B-4    |  |  | \$390.00** |
| <input type="checkbox"/> | SV171/21714RB | Upgrade 17½"x17" Pillow-Top Recline Seat Vinyl |  |  | \$150.00   |
| <input type="checkbox"/> | SC171/21814G  | Upgrade 17½" Deluxe Captain Seat Grey Cloth    |  |  | \$150.00   |
| <input type="checkbox"/> | SV201820RG    | Upgrade 20"x18" Recline Capt Grey Vinyl        |  |  | \$150.00   |
| <input type="checkbox"/> | SP181821RG    | Upgrade Pan seat w/ Capt back recline (18x18)  |  |  | \$150.00   |

| ARM REST ASSY            |              |                              |  |  |            |
|--------------------------|--------------|------------------------------|--|--|------------|
| <input type="checkbox"/> | 104300-88800 | (L) Arm Rest Assy w/ Arm Pad |  |  | \$160.00** |
| <input type="checkbox"/> | 104400-88800 | (R) Arm Rest Assy w/ Arm Pad |  |  | \$160.00** |

Note: Pan seats avail in 16x16 & 20x20 / cushion not sold

| TIRES                    |              |              |                                |    |           |
|--------------------------|--------------|--------------|--------------------------------|----|-----------|
| <input type="checkbox"/> | TAR26085     | Rear         | 10" Pnuematic Tires-all Models | x2 | \$32.00** |
| <input type="checkbox"/> | 205102-88100 | Rear         | 10" Tire Tubes-all Models      | x2 | \$20.00** |
| <input type="checkbox"/> | TAF26085     | Front        | 10" Pnuematic Tires-888B-4     | x2 | \$32.00** |
| <input type="checkbox"/> | 205102-88100 | Front        | 10" Tire Tubes-888B-4          | x2 | \$20.00** |
| <input type="checkbox"/> | TAF280250    | Front        | 9" Pnuematic Tires-888B-3      | x1 | \$36.00** |
| <input type="checkbox"/> | 205202-77100 | Front        | 9" Tire Tubes-888B-3           | x1 | \$20.00** |
| <input type="checkbox"/> | STU8-9(K)    | Front & Rear | 888B-3 Solid Tire Upgrade Kit  |    | \$150.00  |
| <input type="checkbox"/> | STU8A(K)     | Front & Rear | 888B-4 Solid Tire Upgrade Kit  |    | \$200.00  |

| CHARGER                  |              |   |  |  |            |
|--------------------------|--------------|---|--|--|------------|
| <input type="checkbox"/> | 109201-88306 | 3 amp Charger Off-Board (888B-3 & 888B-4) |  |  | \$180.00** |

| ACCESSORIES              |              |                            |  |  |           |
|--------------------------|--------------|----------------------------|--|--|-----------|
| <input type="checkbox"/> | 109901-88101 | Front Wicker Basket        |  |  | \$70.00** |
| <input type="checkbox"/> | CCH-003(K)   | Cane Holder                |  |  | \$150.00  |
| <input type="checkbox"/> | CH-02(K)     | Oxygen Tank Holder (Metal) |  |  | \$150.00  |
| <input type="checkbox"/> | WH-0185(K)   | Walker Holder              |  |  | \$150.00  |
| <input type="checkbox"/> | WIG-0168     | Wig-wag extensions (pair)  |  |  | \$20.00   |
| <input type="checkbox"/> | RB-0125(K)   | Rear Wire Mesh Basket      |  |  | \$170.00  |
| <input type="checkbox"/> | 300909-03    | Safety Belt                |  |  | \$34.00   |

**\*\*STANDARD WITH UNIT**

All specifications are subject to change without notice.  
Shoprider Mobility Products, Inc® reserves the rights of any changes to the unit:

- 1) Include the anti-tip wheel or rear castor.
- 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown below:
  - a) The weight of user.
  - b) Ground surface and terrain.
  - c) Battery capacity and conditions.
  - d) Type of charger.
  - e) Ambient Temperature.
  - f) The way and/or method of driving.
  - g) If the battery and mechanical moving parts are fully broken in.

**TOTAL \$** \_\_\_\_\_

**MIDSIZE SPECIFICATIONS**

**SUNRUNNER-3 (888B-3)**

|                                |                |
|--------------------------------|----------------|
| Base LxWxH*                    | 45"x24"x38"    |
| *Ground to Seat Top            |                |
| Ground Clearance               | 3"             |
| User Weight Capacity           | 300 Lbs        |
| Total Unit Weight              | 165 Lbs        |
| Battery Module Weight (U1)     | 25 Lbs x2      |
| Number, Size Tire (front/rear) | 1,9"/2,10"     |
| BAT-U1 Module Capacity         | 12v35Ah x2     |
| Charger-Off Board              | 3 amp          |
| Turning Radius                 | 37"            |
| Range                          | up to 25 miles |
| Maximum Speed                  | 5 mph          |
| Maximum Climbing Angle         | 14% / 8°       |

**SUNRUNNER-4 (888B-4)**

|                                |                |
|--------------------------------|----------------|
| Base LxWxH*                    | 49"x22"x38"    |
| *Ground to Seat Top            |                |
| Ground Clearance               | 3"             |
| User Weight Capacity           | 300 Lbs        |
| Total Unit Weight              | 180 Lbs        |
| Battery Module Weight (U1)     | 25 Lbs x2      |
| Number, Size Tire (front/rear) | 2,10"/2,10"    |
| BAT-U1 Module Capacity         | 12v35Ah x2     |
| Charger-Off Board              | 3 amp          |
| Turning Radius                 | 50"            |
| Range                          | up to 25 miles |
| Maximum Speed                  | 5 mph          |
| Maximum Climbing Angle         | 14% / 8°       |