

PRICE LIST & ORDER FORM



For easy ordering and information, please contact Customer Service **Toll Free at 1.800.743.0772 FAX: 1.877.797.7081** Effective July 2009



UL8W36 TROOPER

NAME: _____
 PHONE NO: _____
 ACCT NO: _____
 PO #: _____ DATE: _____
 ADDRESS: _____
 CITY: _____ STATE: _____ ZIP: _____
 SHIP TO: _____
 PHONE NO: _____
 ADDRESS: _____
 CITY: _____ STATE: _____ ZIP: _____

SPECIFICATIONS

OVERALL DIMENSIONS

Base Length.....35"
 Base Width.....24"
 Height.....44"
 Ground to Seat Top
 Ground Clearance.....3"

WEIGHTS

User Weight Capacity.....300 Lbs
 Total Weight(w/Batts & Seat).....160 Lbs
 Battery Pack Weight.....55 Lbs
 Captain Seat w/Controller Assy37.1 Lbs

TIRES

Number, Size of Front Tire.....2,6"
 Number, Size of Rear Tire.....2,9"

BATTERIES

BAT-U1 Module Capacity.....12v35Ah x2

CHARGER

Charger-Off Board.....4 amp

MOTOR POWER OUTPUT

Horsepower.....0.45 x2

DRIVING SPECS

Turning Radius32.1"
 Range13-15 mph
 Maximum Speed3.75 mph
 Maximum Climbing Angle10% / 6°

BASE SCOOTER

UL8W36 Trooper Powerchair K0823 \$5,495.00
 Red or Blue

BATTERY

BAT-U1 12V-35AMP x2 \$140.00**ea.

SEATS (W"xD"xH")

SV181/2171/2RG Two-Tone Grey Vinyl Captain Seat \$324.00**
 SP161618RG Upgrade 16"x16" Seat Pan/Captain Back-reclining \$400.00
 Also Available in 18"/ Seat Cushion not included

ARM REST ASSY

104300-00106 (L) Arm Rest Assy w/ Arm Pad \$140.00**
 104400-00106 (R) Arm Rest Assy w/ Arm Pad \$140.00**

TIRES

105201-69400 280X250 Solid Rear Wheel Assy x2 \$112.00**
 1351-00900-00 6" Solid Tire Front x2 \$140.00**

ANTI-TIP WHEELS

209928-89300 Anti-Tip Wheel Assy x2 \$20.00**

CHARGER

109201-88109-A2 3 pin; 4 Amp charger; UL8W; \$230.00**

CONTROLLER

1361-00300-00 Joystick Controller-UL8WPB \$998.00**

FOOT REST

2329-00100-11 Footplate; light weight w. holes \$200.00**
 ELRUL-R(K) Elevating Leg Rest Kit in Red \$232.00
 ELRUL-B(K) Elevating Leg Rest Kit in Black \$232.00

ACCESSORIES

OXY-03 Mesh Oxygen Tank Holder \$60.00
 CCH-003-GK(K) Cane Holder Kit \$150.00
 WH-0185-GK(K) Walker Holder Kit \$150.00

**STANDARD WITH UNIT

TOTAL \$ _____

All specifications are subject to change without notice. Shoprider Mobility Products, Inc® reserves the rights of any changes to the unit:

- 1) Include the anti-tip wheel or rear castor.
- 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown below:
 - a) The weight of user.
 - b) Ground surface and terrain.
 - c) Battery capacity and conditions.
 - d) Type of charger.
 - e) Ambient Temperature.
 - f) The way and/or method of driving.
 - g) If the battery and mechanical moving parts are fully broken in.