



UL8W36

Trooper

PRICE LIST & ORDER FORM Effective November 2006

For easy ordering and information contact Customer Service Toll Free at:

800-743-0772 Fax: 877-797-7081

Name _____	Phone No _____
Acct No _____	PO # _____
Address _____	City/State _____
Ship To _____	Phone No _____
Address _____	City/State _____
	Zip _____

BASE SCOOTER

MSRP

UL8W36 Trooper Powerchair K0823 \$5,495.00

BATTERY

Red or Blue (prejudgement)

BAT-U1 Bat-U1 Sealed Lead Acid x2 \$280.00

SEATS (W"xD"xH")

SV181/2171/2RG Two-Tone Grey Vinyl Captain Seat \$432.00**

SP161618RG Upgrade 16"x16" Seat Pan/Captain Back-reclining \$400.00

Also Available in 18" Seat Cushion not included

ARM REST ASSY

104300-00106 (L) Arm Rest Assy w/ Arm Pad \$140.00**

104400-00106 (R) Arm Rest Assy w/ Arm Pad \$140.00**

TIRES

105201-69400 280X250 Solid Rear Wheel Assy x2 \$112.00**

1351-00900-00 6" Solid Tire Front x2 \$140.00**

ANTI-TIP WHEELS

209928-89300 Anti-Tip Wheel Assy x2 \$20.00**

CHARGER

109201-88109-A2 3 pin; 4 Amp charger; UL8W; \$200.00**

CONTROLLER

1361-00300-01 Joystick Controller-UL8WPB \$998.00**

FOOT REST

2329-00100-11 Footplate; light weight w. holes \$200.00**

ELRUL-R(K) Elevating Leg Rest Kit in Red \$334.00

ELRUL-B(K) Elevating Leg Rest Kit in Black \$334.00

ACCESSORIES

OXY-03 Mesh Oxygen Tank Holder \$60.00

CCH-003-GK(K) Cane Holder Kit \$150.00

WH-0185-GK(K) Walker Holder Kit \$150.00

SMARTIE SPECIFICATIONS

OVERALL DIMENSIONS

Base Length..... 37.0"

w/o footplate..... 25.9"

Base Width..... 23.5"

Height..... 36.5"

(Ground to Seat Top)

Ground Clearance..... 3"

WEIGHTS

User Weight Capacity..... 300Lbs

Total Weight(w/ Bat's & Seat)..... 160 Lbs

Battery pack weight..... 55 Lbs.

Captain seat w/controller assy.... 37.1Lbs

TIRES

Number, Size of Front Tire..... 2, 6"

Number, Size of Rear Tire..... 2, 9"

BATTERIES

BAT-U1 Module Capacity.. (12v35Ah x2)

CHARGER

Charger-Off Board..... 4 amp

MOTOR POWER OUTPUT

Horsepower..... .45hp X2

DRIVING SPECS

Turning Radius..... 32.1"

Range..... 13-15 miles

Maximum Speed..... 3.75 mph

Maximum Climbing Angle..... 10% (6°)

TOTAL \$

**STANDARD WITH UNIT

All specifications are subject to change without notice.

Shoprider Mobility Products, Inc® reserves the rights of any changes to the unit:

- 1) Include the anti-tip wheel or rear castor.
- 2) Driver Weight may exceed weight of the unit; speed must be reduced when turning.
- 3) The actual driving range varies with the factors shown below:
 - a) The weight of user.
 - b) Ground surface and terrain.
 - c) Battery capacity and conditions.
 - d) Type of charger.
 - e) Ambient Temperature.
 - f) The way and/or method of driving.
 - g) If the battery and mechanical moving parts are fully broken in.